



## surgiCam® HD/4K wireless

### Your benefits

**Independence** – can be retrofitted at any time, as no wiring is required

**Simplicity** – easy and quick installation

**Flexibility** – simply take the camera with you into the next room

**Excellence** – professional image quality with maximum sharpness of detail – particularly with the surgiCam® 4K

### Technical data at a glance

Camera Systems	surgiCam® HD wireless	surgiCam® 4K wireless
	<b>High definition</b>	<b>Ultra high definition</b>
Resolution	1920 x 1080	3840 x 2160
Effective pixel number	> 2,000,000	> 8,000,000
Zoom	10 x optical 12 x digital	20 x optical 12 x digital
Motorized image alignment	yes	yes
Video signal	HDMI (DVI with adapter)	HDMI
IP protection class	IP 55	IP 55
Control		
marLED® X	marTouch® Pro, handle, light body	marTouch® Pro, handle, light body
marLED® E	marTouch®, handle, light body	marTouch®, handle, light body
Sterile cover	Single-use cover / reusable	Single-use cover / reusable

## KLS Martin Group

### **KLS Martin Australia Pty Ltd.**

Sydney · Australia  
Tel. +61 2 9439 5316  
australia@klsmartin.com

### **KLS Martin do Brasil Ltda.**

São Paulo · Brazil  
Tel. +55 11 3554 2299  
brazil@klsmartin.com

### **KLS Martin Medical (Shanghai) International Trading Co. Ltd.**

Shanghai · China  
Tel. +86 21 5820 6251  
china@klsmartin.com

### **KLS Martin India Pvt Ltd.**

Chennai · India  
Tel. +91 44 66 442 300  
india@klsmartin.com

### **Martin Italia S.r.l.**

Milan · Italy  
Tel. +39 039 605 67 31  
italia@klsmartin.com

### **Nippon Martin K.K.**

Tokyo · Japan  
Tel. +81 3 3814 1431  
nippon@klsmartin.com

### **KLS Martin SE Asia Sdn. Bhd.**

Penang · Malaysia  
Tel. +604 505 7838  
malaysia@klsmartin.com

### **KLS Martin de México S.A. de C.V.**

Mexico City · Mexico  
Tel. + (52) 55-7572-0944  
mexico@klsmartin.com

### **Martin Nederland/Marned B.V.**

Huizen · The Netherlands  
Tel. +31 35 523 45 38  
infonl@klsmartin.com

### **Gebrüder Martin GmbH & Co. KG**

Moscow · Russia  
Tel. +7 499 792-76-19  
russia@klsmartin.com

### **KLS Martin Taiwan Ltd.**

Taipei 106 · Taiwan  
Tel. +886 2 2325 3169  
taiwan@klsmartin.com

### **Gebrüder Martin GmbH & Co. KG**

Dubai · United Arab Emirates  
Tel. +971 4 454 16 55  
middleeast@klsmartin.com

### **KLS Martin UK Ltd.**

Reading · United Kingdom  
Tel. +44 1189 000 570  
uk@klsmartin.com

### **KLS Martin LP**

Jacksonville · Florida, USA  
Tel. +1 904 641 77 46  
usa@klsmartin.com

### **Gebrüder Martin GmbH & Co. KG**

#### **A company of the KLS Martin Group**

KLS Martin Platz 1 · 78532 Tuttlingen · Germany  
P.O. Box 60 · 78501 Tuttlingen · Germany  
Tel. +49 7461 706-0 · Fax +49 7461 706-193  
info@klsmartin.com · www.klsmartin.com