Arm and Hand
Surgery Tables



Introducing our innovative **arm and hand surgery tables**

CONTENTS

ARM & HAND SURGERY TABLES

THE GARDINER HAND SURGERY SYSTEM™ Seated Position	3
THE MAGNA HAND SURGERY TABLE Supine, Seated Position, Partially Elevated Position	4
THE WOODLAND ARM & HAND SURGERY TABLE Supine	5
ACCESSORIES	6
COMPATIBILITY	7

At Anetic Aid, our mission is to deliver the best outcomes for patients and best experiences for clinical teams. We listen carefully and understand needs.

This allows us to design, develop and manufacture, as well as establishing partnerships, that lead to the supply of innovative, intuitive, high-quality medical devices for patient transport, treatment and recovery.

This booklet allows us to bring you our latest arm & hand surgery tables, The Gardiner, The Magna and The Woodlands. These arm and hand tables are suited for high volume, low complexity hand surgery in numerous patient positions, and are developed in conjunction with leading hand surgeons. We also demonstrate the essential accessories you need and their compatibility with both our QA4™ Mobile Surgery System and our QA3 Stretcher range.

The Gardiner
Hand Surgery System™

The perfect solution for high volume, low complexity hand surgery.

Developed in conjunction with consultant plastic & hand surgeon Matthew Gardiner and Frimley Health NHS Foundation Trust, this hand table is perfect for supporting hand surgery using the WALANT approach - with the patient in a comfortable seated position.

The Gardiner Hand Surgery System™ includes the hand surgery table, flexible screen pole and visibility screen.

Features:

- **Dimensions - L730mm x W335mm**
- **Weight - 2.8kg**
- **Perfect for a huge range of hand surgeries, including de Quervain's release and carpal tunnel decompression**
- **Enables the WALANT approach**



Description	Code
The Gardiner Hand Surgery System™ *	21380

***NOTE: CLAMPS SOLD SEPARATELY**

The Magna
Hand Surgery Table

Supine, seated or partially elevated? The Magna Hand Surgery Table covers it all.

Developed in conjunction with a consultant hand surgeon, this hand table is perfect for day surgery or operating theatre departments where high volume, low complexity hand surgery can be undertaken, following the WALANT approach or when utilising a tourniquet.

Radiolucent, variable height and supplied with a firm surface pad, it can be used when the patient is supine, partially elevated, or in a comfortable seated position.



Features:

- **Dimensions - L740mm x W450mm**
- **Weight - 5.4kg**
- **Covers a huge range of hand surgeries in supine, seated or partially elevated positions**
- **Includes a laterally-sliding, fixed, fold-away telescopic support leg for easier C-arm access.**

Description	Code
The Magna Hand Surgery Table*	21384

***NOTE: CLAMPS SOLD SEPARATELY**

The Woodlands
Arm & Hand Surgery Table

Our specially designed surgery table offering more rigidity in the supine position.

Developed in conjunction with a consultant hand surgeon, The Woodlands is perfect for operating theatre or day surgery departments where arm and hand surgery requires true rigidity for more complex procedures - with the patient supine.

The hand table is radiolucent and is supplied with a firm surface pad.



Features:

- **Dimensions - L690mm x W450mm**
- **Weight - 7.6kg**
- **Offers more rigidity for complex hand surgeries in the supine position**
- **Includes twin integrated direct-on lever lock clamps.**
- **Includes a laterally-sliding, fixed, fold-away telescopic support leg for easier C-arm access.**

Description	Code
The Woodlands Arm & Hand Surgery Table	21385

ACCESSORIES



QA3 Accessory Bracket

A new side rail accessory bracket which can be fitted to our QA3™ Stretchers, allowing for other QA3 accessories to be fitted with ease.

Description	Code
QA3 Accessory Bracket	21156



Clamps

Compatible with our Gardiner™ and Magna Hand Surgery Tables, these clamps provide security and stability for many different procedures.

Description	Code
Clamp - Socket - Stainless Steel - for UK, EU & US Side Bar - **FOR THE MAGNA**	10248
Clamp - Socket - Stainless Steel (10-308-UK) - Hook on with Lever Lock - UK Side Bar - V STUD **- FOR THE GARDINER**	10248-B1
Clamp - Socket - Stainless Steel (10-308-EU) - Hook on with Lever Lock - EU Side Bar - V STUD **- FOR THE GARDINER**	10248-B2

Theatre Pillow with Elasticated Retention Strap

This superbly comfortable, water-resistant theatre pillow comes complete with elasticated strap, helping it retain its position on a QA3™.

Description	Code
Pillow with Elasticated Retention Strap	21188



COMPATIBILITY

Not only compatible with the QA4™ Mobile Surgery System, our new arm and hand surgery tables are also **fully compatible** with our QA3 Stretcher range. See below:



The Gardiner
Hand Surgery System™

- QA4 Mobile Surgery System
- QA3 Patient Stretchers (Patient, DRIVE)
- QA3 Emergency Stretchers (Emergency, DRIVE, Powered)
- QA3 Powered Treatment Stretcher



The Magna
Hand Surgery Table

- QA4 Mobile Surgery System
- QA3 Patient Stretchers (Patient, DRIVE, Treatment)
- QA3 Emergency Stretchers (Emergency, DRIVE, Powered)
- QA3 Powered Treatment Stretcher



The Woodlands
Arm & Hand Surgery Table

- QA4 Mobile Surgery System
- QA3 Patient Stretchers (Patient, DRIVE, Treatment)
- QA3 Emergency Stretchers (Emergency, DRIVE, Powered)
- QA3 Powered Treatment Stretcher

SEE ALSO

AneticAid®



QA4 Mobile Surgery System

The complete multi-purpose solution for transport, treatment and recovery.

PLUS

our other product areas

Scan this code to see our **QA4's** journey through day surgery!









For more information on any of our products or service contracts, call:

sales@aneticaid.com **+44 (0)1943 878647** aneticaid.com

DESIGN & INNOVATION

MADE IN THE UK

GLOBAL DISTRIBUTION

Innovative **Medical** Technology - **Practically Applied**

Anetic Aid Ltd, Sapper Jordan Rossi Park, Otley Road, Baildon BD17 7AX, UK