Introducing the
Split Leg Positioner - with EU Side Bar

Add even greater surgical versatility to your QA4™ with this pair of individually removable, multi-position leg sections that include full length EU sized accessory side bar.

The Split Leg Positioner further enables the use of the QA4™ for a number of surgical procedures within such specialities as gastro-oesophageal, bariatric and orthopaedics - without the need to pre-prepare the QA4™ with them attached, or the need to transfer the patient.

If you want any information on the Split Leg Positioner, QA4™ Mobile Surgery System or QA4™ accessories, please contact your Regional Account Manager or call us on +44 (0)1943 878647.

Want more information?

Scan the code below to discover the QA4™ Mobile Surgery System - the UK's market leading complete multi-purpose solution for patient transport, treatment and recovery.



QA4™

Powered Mobile Surgery System

Split Leg Positioner - with EU Side Bar

For more information on any of our products contact us:

sales@aneticaid.com
+44 (0)1943 878647
aneticaid.com

AneticAid®

AneticAid®



Features & Benefits

Twist grip handles
intuitive and easy to use

Split leg section
offers greater manouvrability and surgical access

Multi-purpose specialities
including gastro-oesophageal bariatric and orthopaedics

No pre-preparation needed
can attach with ease

Greater surgical versatility
for your existing QA4™

Improved infection control
no need to transfer patient

Available accessories
including clamps, foot support and an extremity belt



Products & Accessories

Item No.	Item Description
21359	QA4™ Split Leg Positioner - with EU Side Bar
10250	Clamp - Universal (Socket & Rail) - Handwheel
10460-A9	Foot Support (10-296) - includes Removable Pads
10843	Oasis Extremity Belt (OA113)

Product Specifications

- **Dimensions (overall):** L1000mm x W335mm
- **Dimensions (leg panels):** L650mm x W350mm
- **Abudction range:** +25° / -10°
- **Lithotomy range:** +30° / -50°
- **Weight:** 13kg
- **Maximum patient weight:** 160kg