CARE SURGICAL



CARE SURGICAL BROCHURE - ISSUE 12 v6

CARE SURGICAL

Specialists in Patient Positioning.

Care Surgical is a global company focused on providing safe and comfortable patient positioning products for a variety of surgical procedures. Offering a comprehensive range of disposable and reusable pads, covers and cushions, our products are compatible with most spine and standard operating table systems on the market.

We pride ourselves on being responsive to our customers' needs and welcome new concepts to complement our product range. Our world class customer service and robust quality management systems have earned us ISO 13485 certification.

Our commitment to excellence is why so many customers choose Care Surgical for their patient positioning needs.



Visit www.care-surgical.com or speak to a member of our team for our latest product portfolio.



Contents



CS Universal Spine Extensions™

4-5 CS Universal Spine Extension™
6-7 CS Carbon Spine Extensions™



CS Prone Head Support Systems & Face Cushions

8 CS Prone Head Support Systems
9 Other Prone Head Supports



Chest & Pelvic Supports

10 CS Easy Prone™
11 CS Spine Frame™
12 CS Prone Patient Supports
13 CS Tony Prone™



CS Universal Arm Supports

14 CS Universal Arm Supports™
15 Arm Support Pads



General Table Accessories

16 Table Accessories



Leg & Foot Supports

17 Leg & Foot Supports/ Stirrup Replacement Pads



Lithotomy Stirrups

18-19 CS FlexiFin™ Stirrups



Disposables

20 CS Disposable Covers/ Supports
21 CS Patient Eye Protectors



Supine Supports

22 CS Supine Supports



Operating Table Accessories

23 Operating table accessories

All Care Surgical products are latex free.

The CS Universal Spine Extension™

The CS Universal Spine Extension™ can be fitted to any standard operating table to provide a patient position and imaging capabilities comparable to a specialist dual-rail spine table.

- + **360°** C-Arm/O-Arm access
- + **Improved efficiency** –simply remove the extension to return the table for standard use
- + **Significantly more cost-effective** in comparison to purchasing a full spine table
- + **Powered leg** for ease of adjustment
- + **Light weight** – one person required to set up
- + Adaptor attaches to **side rails**

CS Universal Spine Extension™ Powered Leg



CS Universal Spine Extension™ System



Standard Set Includes:

CS Universal Spine Extension™ - Table Adaptor

Qty: 1 Unit

CS Universal Spine Extension™ - Powered Leg

Qty: 1 Unit

CS Universal Spine Extension™ - Frame

Qty: 1 Unit

CS Universal Spine Extension™ - Set (EU)

Code: CSM2958
Qty: 1 Set

CS Universal Spine Extension ™- Set (US)

Code: CSM3031
Qty: 1 Set

CS Universal Spine Extension ™- Set (AUS)

Code: CSM3035
Qty: 1 Set

CS Universal Spine Extension ™- Set (UK)

Code: CSM3037
Qty: 1 Set

CS Universal Spine Extension ™- Set (Asia)

Code: CSM3033
Qty: 1 Set

Seamless table breaks providing tailored support and flexibility





Leg height -
Max: 960mm/37.7″
Min: 560mm/22″



Store all accessories on the CS Extension Cart

CS Universal Spine Extension™ Cart

Code: CSM2954
Qty: 1 Unit

CS Carbon Spine Extensions™



+ **Improved efficiency** - simply remove extension and replace sections to return to a standard operating table.
+ **Greatly improved anaesthetic access** - no end columns between the anaesthetic team and the patient
+ **250kg/550lbs. weight capacity** – head end, **180kg/400lbs. weight capacity** – foot end



Perfect for surgeries requiring a table break

Table can be flexed and reflexed allowing the intervertebral space to be opened and closed with the hand control.



Lateral Approach made easy

The CS Carbon Fibre Supine Tops allow the patient to be laterally enabling lateral approach surgery such as XLIF



Solution for supine

The CS Carbon Fibre Supine Tops allow the patient to be positioned supine for anterior approach surgery such as ALIF



Pain management support

The CS Carbon Extension can be utilised for pain management along with other procedure which require an imaging table allowing for 360 degree artifact free imaging.



Compatibility

The CS Carbon Spine Extensions are fully compatible with Wilson style frames.



Detachable support leg for added stability

- + **Light weight** – one person required to set up
- + **360 degree C-Arm/0-Arm access**
- + **Significantly more cost-effective** in comparison to purchasing a full spine table

Compatible Tables

Schmitz Diamond Table:

CS Carbon Spine Extension™ Set- Schmitz

Code: CSM2932
Qty: 1 Set

Merivaara Grand Promerix Table:

CS Carbon Spine Extension™ Set- Merivaara

Code: CSM2930
Qty: 1 Set

OPT Vanto II Table:

CS Carbon Spine Extension™ Set- OPT

Code: CSM2931
Qty: 1 Set



CS Carbon Spine Extension™ Cart

Code: CSM3084
Qty: 1 Unit

CS Prone Plus Head Support System™ & Face Cushions

Distributes the pressure evenly across patient's forehead and cheekbones

The chin is free from pressure with both our Prone Head Support System™ and Prone Plus Head Support System™





CS Prone Plus Face Cushion

Code: CSM2512
Qty: 10 Units
(Also available in small: CSM2514)



CS Prone Plus Head Support System™

Code: CSM2565
Qty: 1 System
(Helmet and Mirror)

CS Prone Head Support System™ & Face Cushion

Open chin design allows for rotation of the patient without removing the endotracheal tube



CS Prone Face Cushion

Code: CSM2510
Qty: 10 Units



CS Prone Head Support System™

Code: CSM2560
Qty: 1 System
(Helmet and Mirror)



Other Prone Head Supports



CS ProneOne™ Head Support

Code: CSM2520
Qty: 8 Units



CS ProneOne™ Head Support Small

Code: CSM2519
Qty: 10 Units





CS Prone Face Block 7" – Right Intubation Slot (compression packed)

Code: CSM2535V
Qty: 10 per box



CS Prone Face Block 7"– Left & Right Intubation Slot (compression packed)

Code: CSM2536V
Qty: 10 per box

CS Prone Face Block 5" – Right Intubation Slot (compression packed)

Code: CSM2550V
Qty: 14 per box

For more options to raise the patients head, go to page 15



Allen® is the registered trademark of Allen Medical Systems Inc.



CS Prone C-Flex Compatible Cushion Small

Code: CSM2515
Qty: 10 Units



CS Prone C-Flex Compatible Cushion Large

Code: CSM2516
Qty: 10 Units



CS Prone Face Mask

Code: CSM2530
Qty: 5 Units

CS Easy Prone™

The CS Easy Prone simplifies optimal chest and pelvis positioning with ergonomic components that support spinal alignment, distribute pressure evenly, and protect sensitive areas

Available Sizes

CS Easy Prone™ - 150mm/6"
- Chest & Pelvis Set

Code: CSM2699
Qty: 1 Set

CS Easy Prone™ - 170mm/6.7"
- Chest & Pelvis Set

Code: CSM2918
Qty: 1 Set

CS Easy Prone™ - 200mm/7.8"
- Chest & Pelvis Set

Code: CSM2919
Qty: 1 Set





Three different heights to suit many different body types and head supports

Standard and comfort covers also available





CS Easy Prone™ Comfort Covers

Code: CSM2926
Qty: 6 Sets



CS Easy Prone™ Standard Covers

Code: CSM2929
Qty: 6 Sets

The CS Spine Frame™

This innovative platform fits all standard operating tables, giving you a patient position comparable to that of a specialist spinal table but at a fraction of the cost.



Ideal for use during surgeries requiring a neutral position, for example Laminectomy and Kyphoplasty procedures.

550 lbs. (250kg) patient weight capacity



CS Spine Frame™ Set

Code: CSM2715
Qty: 1 Unit



CS Comfort Covers - 3 Piece

Code: CSM2606-1
Qty: 6 Sets



The design of the supports allows the abdomen to hang freely to reduce venal pressure.



Simply clip the attachments to the frame.



Covers available to protect the supports and aid infection control. 3 sets included.





CS Prone Patient Supports

Code: CSM2716
Qty: 1 - (5 Piece Set)

The CS Prone Patient Supports are made from multiple layers of pressure-relieving foam

+ Velcro-free chest support
+ Quick-release locking mechanism
+ Compatible with all Allen/Jackson - style spine rail tops

Included Accessories



CS Chest Support

Code: CSM2709
Qty: 1 Unit



CS Hip Supports (Pair)

Code: CSM2717
Qty: 1 Pair



CS Thigh Supports (Pair)

Code: CSM2718
Qty: 1 Pair



CS Comfort Covers - 5 Piece

Code: CSM-2600
Qty: 6 Sets (3 Sets Included)

*Compatible with rail sizes that measure 1.25″ (3.2cm) wide and 1.5″ (3.8 cm) tall with 14.5″ (36.8 cm) space between.
Allen® is the registered trademark of Allen Medical Systems Inc.
*Operating table not included in our portfolio.

The CS Tony Prone™



CS Tony Prone™

Code: CSM2906
Qty: 1 Unit

For use with:

+ Pain Management
+ General Surgery
+ Plastic Surgery
+ And other short prone position procedures



Reusable tabletop support for safe and comfortable prone positioning.



Recessed design* allows for improved abdominal decompression.

* Width x height 200mm x 110mm

CS Universal Arm Support

Multi-axis armboard for uncompromised positioning of the patient's arms.
Suitable for use in all positions.







Features

+ One step set-up
+ Central clamping mechanism for easy, fast and secure positioning of the arm support
+ Rapid locking base allows connection to the side rails of the table
+ Uncompromising, efficient handling and safe holding
+ Space saving on the operating table, compact design



CS Universal Arm Support™

Code: CSM2646
Qty: 1 Unit



CS Universal Arm Support - Moulded™

Code: CSM2647
Qty: 1 Unit



CS Universal Arm Support - Extended™

Code: CSM2649
Qty: 1 Unit



CS Armboard Pad 2″

Code: CSM2733-1
Qty: 1 Unit



CS Armboard Pad 3″

Code: CSM2733-2
Qty: 1 Unit



CS Armboard Pad- with Strap

Code: CSM2576-20
Qty: 1 Unit



CS Armboard Pad 2″ - with Strap

Code: CSM2733-4
Qty: 1 Unit



CS Armboard Pad 3″ - with Strap

Code: CSM2733-5
Qty: 1 Unit



CS Armboard Pad

Code: CSM2576-21
Qty: 1 Unit



*For use with CS Universal Arm Support - Extended

CS Armboard L Pad*

Code: CSM2960
Qty: 1 Unit



CS Pad 1″

Code: CSM2568-1
Qty: 1 Unit



CS Pad 2″

Code: CSM2568-2
Qty: 1 Unit



CS Pad 3″

Code: CSM2568-3
Qty: 1 Unit

Rail Clamps



CS Table Rail Adapter

EU Code: CSM2739-1
US Code: CSM2739-3
Qty: 1 Unit



CS Din Rail Clamp

EU Code: CSM2763
US Code: CSM2765
Qty: 1 Unit

Dual Rail Table Accessories



CS Supine Top with 2″ Foam Pad
(12″ x 20″ x 2″)

Code: CSM2722
Qty: 1 Set



CS Supine Top - 3″ Foam Pad
(12″ x 20″ x 3″)

Code: CSM2722-3
Qty: 1 Unit



CS Supine Top - 2″ Foam Pad
(12″ x 20″ x 2″)

Code: CSM2722-2
Qty: 1 Unit



CS Supine Top - Base Only

Code: CSM2722-7
Qty: 1 Unit



CS Head Tray Attachment- Lowered

Code: CSM2751
Qty: 1 Unit



CS Head Tray Attachment

Code: CSM2714
Qty: 1 Unit



Leg & Foot Supports





CS Knee Rings (Pair)

Code: CSM2821
Qty: 1 Pair



CS Easy Shin Support™

Code: CSM3020
Qty: 1 Unit

CS Replacement Stirrup Pads

+ High-quality, wipe-clean material for durability and improved infection control
+ 20mm/0.7″ memory foam for improved patient comfort and support





CS Yellofin Compatible Replacement Pads

Code: CSM2734
Qty: 1 Pair



CS Ultrafin Compatible Replacement Pads

Code: CSM2741
Qty: 1 Pair



CS Standard Replacement Stirrup Pads

Code: CSM2750
Qty: 1 Pair

CS FlexiFin™ & FlexiFin™ Max Stirrups



CS FlexiFin™ Max Stirrups (Pair) - 360kg/790lbs

Code: CSM2913
Qty: 1 Pair

CS FlexiFin™ Stirrups (Pair) - 227kg/500lbs

Code: CSM2754
Qty: 1 Pair

The CS FlexiFin™ Stirrups are designed for effortless, optimal lithotomy positioning.

+ Flexible Boot to accommodate larger patient calves
+ Pressure-relieving padding reduces risk of pressure injury and compartment syndrome
+ Longer handle for ease of adjustment
+ 227kg (500lbs)/360kg (790lbs) max patient weight, so suitable for a wide range of patients









CS FlexiFin™ & FlexiFin™ Max Stirrups - Hook & Loop

CS FlexiFin™ Max Stirrups - Hook & Loop (Pair)

Code: CSM2756
Qty: 1 Pair

CS FlexiFin™ Stirrups - Hook & Loop (Pair)

Code: CSM2755
Qty: 1 Pair

No fin for Improved access



CS FlexiFinless™ Stirrups

CS FlexiFinless™ Stirrups (Pair)

Code: CSM3050
Qty: 1 Pair

CS FlexiFinless™ Stirrups - Hook & Loop (Pair)

Code: CSM3058
Qty: 1 Pair

Stirrup Accessories



CS Stirrup Cart

Code: CSM2724
Qty: 1 Unit



CS Blade Clamp – UK/EU

Code: CSM2761
Qty: 1 Unit



CS Blade Clamp Plus

US Code: CSM2762
EU Code: CSM2767
Qty: 1 Unit

Features

Patient Weight Capacity	Standard: 227kg/500lbs. Max: 360kg/790lbs.
Lithotomy Adjustments (From sterile field)	Yes
Lithotomy Range (From sterile field)	+90° to -45°
Abduction Range (From sterile field)	+25° to -10°
Maximum medial/ lateral shift/boot rotational angle	+-23°

Comfort/Standard Covers



CS Comfort Covers- 3 Piece

Code: CSM2606-1
Qty: 6 Sets



CS Comfort Covers - 5 Piece

Code: CSM2600
Qty: 6 Sets



CS Comfort Covers Large - 5 Piece

Code: CSM2603
Qty: 6 Sets



CS Pro-Axis Compatible Set

Code: CSM4570
Qty: 6 Sets



CS Wingset Comfort Covers

Code: CSM-2605
Qty: 5 Sets



CS Bow & Wilson Comfort Covers

Code: CSM2608
Qty: 6 Sets



CS Bow & Wilson Standard Covers

Code: CSM2570
Qty: 5 Sets



CS Jackson Comfort Covers - 5 Piece

Code: CSM2602
Qty: 6 Sets



CS Jackson Standard Covers - 5 Piece

Code: CSM2610
Qty: 6 Sets

Other disposables



CS Disposable Leg Holder

Code: CSM2738
Qty:12 Units



CS Head Strap Disposable

Code: CSM2667
Qty: 10 Units



CS Traction Boot Foam Insert

Code: CSM2666
Qty: 6 Units



CS Disposable Lateral Head Support (Adult)

Code: CSM2410
Qty: 20 Units



CS Memory Foam Sheet

Code: CSM2588
Qty: 2 Units



CS Egg Box Foam 80" x 20"

Code: CSM2585
Qty: 1 Unit



CS Standard Arm Cushions

Code: CSM2575
Qty: 10 Pairs



CS Torso Support Cushion - Medium

Code: CSD-150-M
Qty: 2 Units



CS Adjustable Prone Leg Cushion

Code: CSM2590
Qty: 2 Units



CS Torso Support Cushion - Small

Code: CSD-150-S
Qty: 2 Units



CS Prone Arm Secure - Large

Code: CSM2674
Qty: 10 Units

CS Patient Eye Protectors

The soft foam CS Patient Eye Protector protects the delicate nerves and skin around the patient's eyes and can be adapted to every face.



+ Single-use to aid infection control
+ Reduces risk of pressure-related injury
+ Protects the eyes from damage and dehydration



CS Patient Eye Protectors

Code: CSM2790
Qty: 25 Units

CS Plexus Wedges



CS Plexus Wedges (Pair)

Code: CSM2908
Qty: 1 Pair



Placed under the patients shoulder blades in supine position the wedges help to reduce the pressure on the brachial plexus and to improve the comfort e.g. different degrees of a hunched back.

CS Supine Bariatric Support

The Supine Bariatric Support is a 2-piece wedge cushion set to elevate the head and thorax in the supine position, reducing pressure and facilitating intubation.

+ Helps to easily reach the optimal intubation position by raising the head, upper body and shoulders above the chest (“ramped position”)
+ Fully sealed for a better cleaning
+ Fully radiolucent
+ Maximum Weight: 250kg/550lbs.



CS Supine Bariatric Support

Code: CSM2648
Qty: 1 Set

Other Supine Supports



CS Plexus Wedges (Pair)

Code: CSM2908
Qty: 1 Pair



CS Universal Supine Head Support

Code: CSM3024
Qty: 1 Unit



CS Disposable Supine Head Support


Code: CSM2411
Qty: 20 Units



CS Supine & Lateral Head Support

Code: CSM2643
Qty: 1 Unit

*Also available in paediatric size: CSM2644



Operating Table Accessories

Clamps / Rails



Round Bar Accessory Clamp

Code: 10-301
Qty: 1 Unit



Accessory clamp, Non Rotatable

Code: 10-302
Qty: 1 Unit



Radial Clamp Universal

Code: 10-308
Qty: 1 Unit

Lateral Supports



Lateral Support narrow

Code: 10-566
Qty: 1 Units



Rotatable Lateral support narrow

Code: 10-560
Qty: 1 Unit



Rotatable Lateral support wide

Code: 10-561
Qty: 1 Unit



Pubis Support

Code: 10-567
Qty: 1 Unit



Supporting Cushion Vertical

Code: 10-368
Qty: 1 Units

Arm and Hand Rests



Arm & Hand Surgery Table

Code: 10-390
Qty: 1 Unit



Arm & Hand Surgery Table Carbon Fiber

Code: 10-396
Qty: 1 Unit



Arm & Hand Surgery Table Carbon Fiber short

Code: 10-396-K
Qty: 1 Unit



Arm & Hand Surgery Table short

Code: 10-390-K
Qty: 1 Units

Straps



Leg & Body Strap wide with velcro

Code: 10-464-M
Qty: 1 Unit

Anesthesia Screen



Anesthesia Screen for Operating Table

Code: 10-310
Qty: 1 Unit



Anesthesia Frame flexible with clamps

Code: 10-311
Qty: 1 Unit

View the full range of operating table accessories:

Notes

Notes

Notes

Notes

CARE SURGICAL

UK: Care Surgical Ltd, Unit 6 Ringtail Road
Burscough, Lancashire L40 8JY UK
Tel: +44 (0) 1704 336671
info@care-surgical.com
www.care-surgical.com

US: Care Surgical LLC,
43 Broad Street, Suite B106,
Hudson, MA 01749, USA
Tel: 866-243-4107
sales@care-surgical.com
www.care-surgical.com

