




Dr. Mach GmbH & Co. KG
Am Brucker Feld 4,
D-85567 Grafing
Telephone: +49 (0) 8092 / 2093-0
e-mail: info@dr-mach.de
www.dr-mach.de

DR.MACH

Mobile and stable.

Portable and power independent – the new Dr. Mach battery stand



Consistent light supply

The new battery-powered stand takes mobility to a new level – thanks to perfect weight distribution and a high-performance battery, even when there is inconsistent electricity supply. Available in most Dr. Mach LED 130 and LED 150 variations. This offers the perfect solution for mobile medical lighting: portable, stable and cordless.

Maximum flexibility


The Dr. Mach battery stand offers an ideal solution when maximum flexibility is required. Thanks to a special base construction with four smooth-running castors, it is portable and offers a stable stand at the same time. Thanks to the new high-performance battery, which is fixed to the base of the stand, medical lighting is guaranteed even without an external power source.

Technical Data*	
Input voltage	100 - 240 V AC
Output voltage	24 V DC
Frequency	50/60 Hz
Charging time	2,5 hours
Battery running time (depending on model)	approx. 4-8 hours
Interchangeable battery	Yes

*technical data subject to changes

Wide variety of applications

For use within a variety of rooms, without a power source and for luminaire mobility without installation: the Dr. Mach battery stand is the ideal solution for applications where a fixed installation is not possible. Thanks to its specially designed base, precise suspension arm technology and modern rechargeable battery technology, the battery stand offers optimum functionality. It can be used reliably in different rooms, in field hospitals or mobile clinics. The stand can be folded and disassembled for easy transportation. This allows it be stowed away compactly and easily (e.g. in Zarges B20 Y boxes).



Modern Technology

Dr. Mach medical lights feature the latest technologies – for over 75 years. Specializing in the development, construction and use of medical lights, Dr. Mach excels with their reliability and high functionality. Doctors and medical professionals around the world work with lights such as the LED 130 and LED 150, both of which are available in different versions – ideal for various medical facilities. All Dr. Mach products are subject to the highest quality requirements and are Made in Germany.

