KLS martin GROUP



independant® Neo

Flexibility in a new design

Modular ceiling supply units for
operating theaters and intensive care units

independant® Neo – the systems at a glance

Our newly designed independant® Neo console combines maximum flexibility with adaptability and is specifically designed to meet the requirements of modern operating theaters and intensive care units. It offers a future-oriented solution that can be seamlessly adapted to your individual needs while meeting the highest standards in medical technology.


It can be customized to meet your needs. Whether gas, electric or data sockets – choose the equipment that best suits your requirements and enhance it with additional accessories.

In combination with the proven support arms, the independant® Neo console is an efficient and reliable supply solution that has been specifically developed for use in operating theaters and intensive care units.

The systems are available in different versions – with standard support arms for horizontal positioning or with additional height adjustment by means of a motor or spring arm.


The LED BrakeGuide ensures simple and safe handling for systems with electromagnetic breaks. As an alternative option, the system with pneumatic brakes offers equally reliable and precise control. Here, the brakes are identified by color coding.






Independant® Neo S

The horizontally pivoting support arm of the independant® Neo S is designed for a load capacity of up to 1000 kg. This ensures a high degree of flexibility and functionality for the entire system.



independant® Neo M

With its integrated motor, this model offers height adjustment in addition to horizontal positioning. It features an impressive load capacity of up to 250 kg and is perfect for creating efficient and safe working environments.




independant® Neo F

You can also rely on our model with a height-adjustable spring arm: With a load capacity of up to 180 kg, the system offers flexible and reliable support in demanding work environments.

All information depends on the selected configuration.


independant® Neo – the vertical supply console

In a cramped medical environment, every inch counts. The vertical design of the independant® Neo console maximizes usable space by securely and ergonomically accommodating devices, monitors and other essential components in the smallest of footprints. This allows workflows to be made more efficient.









Advantages and features of the vertical independant® Neo console

- **Flextrack system**
  The standard-fitted Flextrack system allows for the stepless, flexible mounting of accessories, rails or handles – entirely according to your needs. This allows accessories to be added at any time. As a result, you remain flexible and prepared for future requirements.

- **Indirect light (optional)**
  Integrated lighting in the base of the console, on the support arms and on the side of the console ensures optimal visibility and a pleasant working atmosphere.

- **Integrated cable management (optional)**
  A seamless cable management system enables all cables to be routed cleanly and securely along the entire length of the console, and can be mounted on both sides if desired.

- **Multifunction rack (optional)**
  Available in different lengths, the multifunction rack offers flexible options for mounting trays and other accessories. With its high load capacity of up to 180 kg, the system offers sufficient capacity for a variety of devices and accessories.

- **Mono-multifunction rack (optional)**
  Option to install a mono MFR system for specific applications. This single MFR can be installed on the opposite side.

- **MediSound system (optional)**
  The only medically approved audio system integrated into ceiling supply units for use in operating and intensive care areas that ensures precise and secure audio transmission via Bluetooth.








All information depends on the selected configuration.

independant® Neo – the new horizontal supply console

The horizontal design of the independant® Neo console also offers medical facilities an excellent opportunity to make their workflows more efficient. With the flexible option of installing outlets on five different sides, this console optimally adapts to the specific requirements of everyday clinical routines.

The major advantage of the horizontal console becomes apparent when combined with height-adjustable support arms. These allow valuable space to be created in cramped operating theaters. Depending on the room height, the console can be moved into a parking position that allows walking underneath the console. This ensures a freer working environment and easier access to important areas – especially in critical situations where every move counts.






Advantages and features of the horizontal independant® Neo console

- **Operation via handle or remote control**
  Whether operated manually using a handle or conveniently by remote control - the independant® Neo horizontal console offers intuitive control that can be easily integrated into daily workflows.

- **Multifunction rack can be retrofitted in different lengths (optional)**
  To provide even greater versatility, the console can be expanded with a retrofittable multifunction rack (load capacity up to 165 kg) in different lengths. This option allows for optimal organization of devices and accessories, tailored to the specific needs of your facility. Allowance for this must be taken into account when placing the order.

- **Mounting of up to three trays (optional)**
  The option of mounting up to three trays provides additional storage space and work surfaces, keeping important instruments and materials within easy reach at all times and promoting a more efficient workflow.

All information depends on the selected configuration.



independant® Neo – numerous combination options

Our independant® Neo ceiling supply units offer you a selection of different combination options that are perfectly tailored to the specific requirements of the different medical fields.

In addition to the classic single and tandem installations of the systems, you can also benefit from a tandem installation that can be made even more flexible by adding a KLS Martin monitor arm or a operating light, saving you the need for additional mounting points.

Additionally, it is possible to install a marLED® E3 examination light on the upper support arm of the ceiling supply unit. This option is particularly useful in intensive care, where it ensures precise illumination during examinations and treatments.

Our solutions are designed to optimize and thus facilitate the workflow of medical personnel.

Opt for flexibility and quality – with our individually adaptable ceiling supply units, which are precisely tailored to your needs.




All information depends on the selected configuration.

independant® Neo F and CEMOR® – flexible monitor suspension systems

In combination with our high-quality support arms, our CEMOR® monitor suspension system offers the perfect solution for the installation of multiple or particularly large monitors in the operating theater. Owing to its robust and flexible design, the system can be optimally integrated into any OR environment, supports precise visualization and contributes to a smooth workflow.

It is available in different versions: as a 4-, 6- or 8-fold variant, as well as being specifically designed for large screens. This versatility allows you to choose the ideal configuration for your specific requirements, thus further increasing the efficiency and quality of your work.



Number of monitors	Possible sizes of monitors (depending on weight and height)	Monitors horizontal	Monitors vertical
4-fold	19“ - 32“	2	2
6-fold	19“ - 32“	3	2
8-fold	19“ - 24“	4	2
LS (Large Screen)	37“ - 60“	1	1

independant® Neo – supplementary documents and products

For more information on the available systems and compatible accessories, as well as complementary products from our portfolio, please refer to the documents below. If you have any questions or require individual advice, please do not hesitate to contact your Area Sales Manager, who will be pleased to assist you at all times.







All information depends on the selected configuration.

KLS Martin Group

KLS Martin Australia Pty Ltd
Sydney · Australia
Tel. +61 2 9439 5316
australia@klsmartin.com

KLS Martin do Brasil Ltda.
São Paulo · Brazil
Tel. +55 11 3554 2299
brazil@klsmartin.com

KLS Martin Medical (Shanghai) International Trading Co., Ltd
Shanghai · China
Tel. +86 21 5820 6251
info@klsmartin.com

KLS Martin India Pvt Ltd.
Chennai · India
Tel. +91 44 66 442 300
india@klsmartin.com

KLS Martin Italia S.r.l.
Milan · Italy
Tel. +39 039 605 67 31
info@klsmartin.com

KLS Martin Japan K.K.
Tokyo · Japan
Tel. +81 3 3814 1431
japan@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.
Penang · Malaysia
Tel. +604 261 7060
malaysia@klsmartin.com

KLS Martin de México, S.A. de C.V.
Mexico City · Mexico
Tel. +52 55 7572 0944
mexico@klsmartin.com

KLS Martin Nederland B.V.
Huizen · Netherlands
Tel. +31 35 523 45 38
infonl@klsmartin.com

KLS Martin Poland Spółka z o.o.
Wrocław · Poland
Tel. +48 605 726 069
info@klsmartin.com

KLS Martin SE & Co. KG
Moscow · Russia
Tel. +7 499 792 76 19
russia@klsmartin.com

KLS Martin Taiwan Ltd.
Taipei · Taiwan
Tel. +886 2 2325 3169
taiwan@klsmartin.com

KLS Martin SE & Co. KG
Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

KLS Martin UK Ltd.
Reading · United Kingdom
Tel. +44 118 467 1500
info.uk@klsmartin.com

KLS Martin LP
Jacksonville · Florida, USA
Tel. +1 904 641 77 46
usa@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.
Hanoi · Vietnam
Tel. +49 7461 706-0
vietnam@klsmartin.com

KLS Martin SE & Co. KG
A company of the KLS Martin Group
KLS Martin Platz 1 · 78532 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com

91-350-62-02 · Rev. 01 · 2026-08 · Printed in Germany · Copyright by KLS Martin SE & Co. KG · All rights reserved · Subject to technical changes