KYRA MEDICAL

Your trusted partners in patient positioning

PATIENT POSITIONING PORTFOLIO



2026
PRODUCT CATALOG

Who We Are

KYRA Medical creates innovative patient positioning solutions for the operating room, designed to minimize risk and optimize patient positioning for improved surgical outcomes across all specialties including gynecology, urology, laparoscopy, colorectal, general, spine, and robotic surgery.

Already known for our iconic purple stirrups, KYRA's design DNA focuses on products that are safer, easier to use, and better for both patients and clinicians. Our range of high-quality, affordable solutions aims to simplify and standardize optimal patient positioning across the surgical suite.

KYRA Quick Facts

2	5	30+	19+	50+	50+
Global HQs	International Service Centers	Patents Held	US Distributors	Global Distributors	OEM Partners



MADE IN THE USA

KYRA products are designed and manufactured in the USA and UK. With dual headquarters and regional service centers in France and Japan, customers are supplied from the optimal location, ensuring reliable delivery, seamless support, and convenient access to local servicing and repairs.

Offices in Boston, USA & Preston, UK


Service & Repair Centers in Boston, USA; Preston, UK; Forbach, France; Kobe, Japan; and Tokyo, Japan.

WHY WE EXIST

To Build Something Better™

KYRA unites a team with decades of combined experience in patient positioning, including those behind the first-ever lithotomy stirrups and innovative spine platforms brought to market.

We prioritize long-term partnerships with clinicians, hospitals, educational institutions, suppliers, and over 30 global medical device leaders. These collaborations inform practical, innovative solutions that improve workflow, enhance safety, and support better surgical outcomes.



The KYRA Design DNA

Powered by Innovation

KYRA Innovation extends beyond our portfolio—our legacy is woven throughout the global medical device industry. Our extensive intellectual property spans proprietary technologies and OEM contributions developed for or licensed by the world's leading medical device manufacturers.

By combining design expertise with real-world clinical insight, we create solutions that are functional, intuitive, and adaptable—advancing patient positioning across our own portfolio and throughout the wider medical device industry.

↑ View our patents

For Patient and Clinician.

Designed for both form and function, KYRA develops innovative, intuitive solutions that support optimal surgical outcomes; engineered to be safer for both patient and clinicians, and thoughtfully designed to be pleasing to the hand and the eye.

Easy-Squeeze Handle
Designed with a reduced squeeze angle and grip force for clinician comfort

Secure-Lok™ Silicone Straps
Durable, secure and easy-to-clean for improved infection control



Purple Touchpoints
Intuitive to use and adjust

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

About Us

Product & Ordering Information

PHONE
U.S. Toll-Free 888-611-KYRA (5972)
Office +1-508-936-3550

FAX
+1-617-977-KYRA (5972)

EMAIL
customers@kyramedical.com

WEBSITE
kyramedical.com

ADDRESS
KYRA Medical, Inc.
102 Otis Street
Northborough, MA 01532 USA

SCAN to visit our website

KYRA COMPLETE SERVICING & WARRANTY

KYRA backs its quality products with superior service and warranty programs designed to keep your operating room running smoothly and confidently.

Standard Warranty
- Hardware: 2 year
- Pads: 1 year

Repair on-site
with specialty kits

Extended Warranty
Product dependent

See full details on KYRAComplete™ p. 80 →

Scan or email for repairs and servicing:
service@kyramedical.com

Brand and Trademark Information

KYRA & KYRA Medical are registered trademarks of KYRA Medical, Inc.

KYRA Comfort, KYRA Clean, KYRA Classic, KYRA Curve, Secure-Lok, Connect, Blade-Lok, Secure Release, ProneAlone, KYRAComplete, “The Purple Stirrup”, “from the creators of the purple stirrup”, “Build something better”, “Your trusted partners in patient positioning", are trademarks of KYRA Medical, Inc. or its subsidiaries.

Allen® and C-Flex™ are registered trademarks of Baxter.

PRODUCT STANDARDS

MADE IN USA

KYRA products are designed and manufactured near Boston, MA USA & Preston, UK.

LATEX

KYRA products are not made with natural rubber latex.

Products may be covered by one or more issued or pending patents.

ISO 13485:2016

EU-MDR

NSF-ISR Registered to ISO 13485

Ordering Information Guide

Each product clearly identifies included components, items required for use, and related or optional accessories that may be added based on clinical needs.

WHAT'S INCLUDED (Example)



Comfort 500 Stirrups
w/Secure-Lok™ silicone straps
KYRA5100

OR



Comfort 500 Stirrups
w/hook and loop straps
KYRA5000

WHAT'S REQUIRED



Secure Release™ Clamp*, p. 40
*2 required
KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

WHAT'S OPTIONAL



Purple™ Stirrup Cart, p. 37
Safely stores and transports stirrups.
KYRA1000



Clean™ Catcher Drainage System, p. 34
All-in-one, disposable urology drainage system.
KYRA8452 (Non-sterile)
KYRA8451 (Sterile)



Urology Bands & Drainage Bags, p. 36
KYRA8400 (Band)
KYRA8410 (Non-sterile)
KYRA8411 (Sterile)

Surgical Position Guide

Icons throughout the catalog indicate compatible surgical positions.



UNIVERSAL



LITHOTOMY



TRENDELENBURG



REVERSE TRENDELENBURG



SUPINE



KIDNEY



LATERAL



PRONE



WILSON




JACKKNIFE

Table Compatibility Guide

Standard OR Table Rails: Alternate order numbers are provided for country-specific rail versions.

Rail Type	Height		Width	
AMERICAN (US)	1.125" (1⅛")	28.6 mm	0.375" (⅜")	9.5 mm
ESCHMANN (UK)	1.25" (1¼")	31.75 mm	0.25 (¼")	6.35 mm
EUROPEAN (EU)	0.98"	25.0 mm	0.40" (~⅜")	10.0 mm
JAPANESE (JP)	1.25" (1¼")	32.0 mm	0.35" (⅜")	9.0 mm
DENYERS (DY)	1.5" (1½")	38.1 mm	0.25 (¼")	6.35 mm

Specify Rail Type (Example)
US KYRA2000 UK KYRA2006
EU KYRA2005 DY KYRA2007
JP KYRA2008



DUAL RAIL COMPATIBLE

Dual Rails: Dual-Rail icon indicates compatibility with spinal operating room tables with 1.25" (3.2 cm) wide and 1.5" (3.8 cm) tall rail sizes.

Product & Ordering Information

Table of Contents

Stirrups/Leg Holders

Stirrup Overview 7
Classic™ Stirrups..... 8
Comfort™ Stirrups..... 18
Clean™ Stirrups 26
Specialty Stirrups 32
Stirrup Storage 37
Leg Holders 38

Urology

Clean™ Catcher Drainage System 34
Urology Bands & Drainage Bags 36

Clamps

Flat Blade Clamps 40
Round Post Clamps..... 41
Side Rail Extenders..... 41

General Accessories

Table Accessories..... 42
Armboards & Pads 44
Arm Positioners 46
Orthopedics..... 48
Patient Transfer 49
Safety Straps & Supports..... 50
Lateral Positioning..... 52

Spine

Curve™ Spine Frame..... 54
Connect™ Spine System 56
Connect™ Head Supports 62
Connect™ Torso Supports..... 66
Connect™ Arm Supports..... 68
Connect™ Supine & Lateral 70
Connect™ Clamps..... 72
Connect™ Accessories 73

Disposables

Head Supports..... 74
Body..... 76

Warranty & Indexes

KYRAComplete™ Warranty 80
Full Product Index..... 82

FIND YOUR FIT

KYRA STIRRUPS

KYRA Classic™

A traditional squeeze-handle stirrup with next-gen innovation.


p. 8



KYRA Comfort™

Widest range of lithotomy for robotics and low lithotomy procedures.

p. 18



KYRA Clean™

Fully fabric-free stirrup with wipe-clean, molded foam pads designed for maximum infection control.


p. 26



KYRA Specialty Stirrups & Accessories

Non-finned boot options for smaller patients.

p. 36



Urology Drainage Solutions

Fully disposable and reusable solutions for fluid drainage.

p. 34




Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

KYRA® Classic™ Stirrups

Reinventing the Classic

Combining patented KYRA technology with classic stirrup design—the Classic stirrup features an improved easy-squeeze handle and traditional straight rod, and can treat your tallest to shortest patients.



Finned and Non-Finned
boot styles available.

Fabric-free silicone straps
help improve infection control.

Straight Rod
Enhanced straight rod classic design.

Easy-squeeze Handle
Reduced squeeze angle and grip force suits smaller hands.

Unique Design
Supports safe positioning of the tallest to shortest patients.



Less grip force needed

Looking for a **fabric-free** stirrup?
Learn about **KYRA Clean™ Stirrups**, p. 26

Premium Servicing & Warranty

KYRA COMPLETE

Replacement Pads

Refresh your stirrups without replacing the entire system. When stirrups remain sound, new boot pads restore comfort and protection.

KYRAComplete™ Extended Warranty

Cover your hardware and prevent downtime with KYRAComplete premium warranty and service packages. KYRA Stirrups are designed to support common repairs on-site with easy repair kits. Details, p. 80.

KYRA0100

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.


Classic Features meet KYRA innovation




✓ Reduced squeeze angle and grip force suits smaller hands



✓ Metal construction supports smooth free-floating boot for patient safety




✓ Staple-free design for a safer surgical environment



✓ Easy-clean, fabric-free silicone SecureLok™ straps help improve infection control

On-site repairs

✓ Maximize uptime with repair kits, extended warranty and premium servicing



✓ Fits a wider range of patient heights and anatomy

Find Your Fit

Model Comparison Guide



	Classic 350	Classic 500	Classic 500 Non-Fin	Classic 500 Plus	Classic 800 Plus
BOOT SIZE/STYLE	Standard Boot with Lateral Fins		Non-finned Std. Boot	Large Boot with Lateral Fins	
PATIENT WEIGHT CAPACITY	350 lbs. (159 kg)	500 lbs. (227 kg)	500 lbs. (227 kg)	500 lbs. (227 kg)	800 lbs. (363 kg)
STIRRUP WEIGHT	13.2 lbs. (6 kg)	13.2 lbs. (6 kg)	14.6 lbs. (6.6 kg)	14.6 lbs. (6.6 kg)	14.6 lbs. (6.6 kg)
BOOT PADS THICKNESS	¾" (19 mm) pads	¾" (19 mm) pads	¾" (19 mm) pads	1" (25 mm) pads	1" (25 mm) pads
STRAP TYPE	Secure-Lok*	Secure-Lok*	Secure-Lok*	Secure-Lok*	Secure-Lok*
LITHOTOMY RANGE	+87° to -35°	+87° to -35°	+87° to -35°	+87° to -35°	+87° to -35°
USE CASE	General lithotomy.		Non-finned boot for robotics access.	General lithotomy— bigger boot size for larger patients.	
PRODUCT NO.	KYRA3800	KYRA5800	KYRA7800	KYRA5850	KYRA8850

All Classic Stirrups are sold in pairs with pads included. Requires two Secure Release™ Clamps.
***Secure-Lok Straps are made of durable silicone, eliminating use of hook and loop.**

Complete ordering details →

KYRA Classic™ 350 Stirrups

FOR PATIENTS UP TO 350 LBS. (159 KG.)

The Classic stirrup features an improved easy-squeeze handle and straight rod. The patented design allows maximum rod length to accommodate the widest range of patient heights. Learn about the full line, p. 8.



Fabric-free silicone straps
help improve infection control.

Easy-squeeze Handle
Reduced squeeze angle and grip force suits smaller hands.

Unique Design
Patented design supports safe positioning of the tallest to shortest patients.

Shown with KYRA Clean™ Catcher Drainage System, p. 34

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Classic™ 350 Stirrups

KYRA3800

0.75" (19mm) thick pads included. Sold in pairs.
Patient Weight Capacity: 350 lbs. (159 kg)
Stirrup Weight: 13.2 lbs. (5.9 kg)
Lithotomy Range: +87° to -35°
Abduction Range: +25° to -9°
Boot Rotation Angle: +/- 25°

- Easy-squeeze handle with reduced squeeze angle and grip force suits smaller hands.
- Position your shortest to tallest patient.
- Easy-clean silicone straps designed for improved infection control.
- Larger boot with ¾" (19 mm) of pressure-reducing padding for modern patient anatomy.



Fabric-free silicone straps

Premium stainless steel gas shock

WHAT'S INCLUDED



Classic 350 Stirrups
w/Secure-Lok™ silicone straps
KYRA3800

WHAT'S REQUIRED



Secure Release™ Clamp*, p. 40
*2 required
KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

WHAT'S OPTIONAL



Purple™ Stirrup Cart, p. 37
Safely stores and transports stirrups.
KYRA1000



Clean™ Catcher Drainage System, p. 34
All-in-one, disposable urology drainage system.
KYRA8452 (Non-sterile)
KYRA8451 (Sterile)



Urology Bands & Drainage Bags, p. 36
KYRA8400 (Band)
KYRA8410 (Non-sterile)
KYRA8411 (Sterile)

POSITION GUIDE

For larger weight capacity, see the
KYRA Classic 500 stirrups →

KYRA COMPLETE Premium Servicing & Warranty

Classic 350 Replacement Pads
Extend the lifetime of your stirrup hardware with replacement stirrup pads.
KYRA6350-P (Secure-Lok Silicone)

KYRAComplete™ Extended Warranty
Cover your hardware and prevent downtime with KYRAComplete premium warranty and service packages. KYRA Stirrups are designed to support common repairs on-site with easy repair kits. Details, p. 80.
KYRA0100

KYRA Classic™ 500/500 Plus Stirrups

FOR PATIENTS UP TO 500 LBS. (227 KG.)

Available with standard boot (Classic 500) or large boot (Classic 500 Plus) to suit a wide range of patient anatomies. The Classic stirrup features an improved easy-squeeze handle and straight rod. Expertly designed to accommodate your tallest to shortest patients. Learn about the full line, p. 8.



Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Classic™ 500 Stirrups

KYRA5800 - Classic 500
KYRA5850 - Classic 500 Plus

Pads included. 500: 0.75" (19 mm); 500 Plus: 1" (25 mm)
Sold in pairs.
Patient Weight Capacity: 500 lbs. (227 kg)
Stirrup Weight: 500: 13.2 lbs. (5.9 kg) • 500 Plus: 14.6 lbs. (6.6 kg)
Lithotomy Range: +87° to -35°
Abduction Range: +25° to -9°
Boot Rotation Angle: +/- 25°



- Easy-squeeze handle with reduced squeeze angle and grip force suits smaller hands.
- Position your shortest to tallest patient.
- Easy-clean silicone straps designed for improved infection control.
- Larger boot sizes with pressure-reducing padding for modern patient anatomy.
- 500 Plus boot option accommodates your largest patients.

WHAT'S INCLUDED



Classic 500 Stirrups
w/Secure-Lok™ silicone straps
KYRA5800

OR



Classic 500 Plus Stirrups
w/Secure-Lok™ silicone straps
KYRA5850

WHAT'S REQUIRED



Secure Release™ Clamp*, p. 40
*2 required

KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

WHAT'S OPTIONAL



Purple™ Stirrup Cart, p. 37
Safely stores and transports stirrups.
KYRA1000



Clean™ Catcher Drainage System, p. 34
All-in-one, disposable urology drainage system.
KYRA8452 (Non-sterile)
KYRA8451 (Sterile)



Urology Bands & Drainage Bags, p. 36
KYRA8400 (Band)
KYRA8410 (Non-sterile)
KYRA8411 (Sterile)

POSITION GUIDE

For a non-finned boot, see the KYRA Classic 500 Non-Fin Stirrup →

KYRA COMPLETE Premium Servicing & Warranty


Classic 500/500 Plus Replacement Pads
Extend the lifetime of your stirrup hardware with replacement stirrup pads.
500: KYRA6470-P (Secure-Lok Silicone)
500 Plus: KYRA6360-P (Secure-Lok Silicone)

KYRAComplete™ Extended Warranty
Cover your hardware and prevent downtime with KYRAComplete premium warranty and service packages. KYRA Stirrups are designed to support common repairs on-site with easy repair kits. Details, p. 80.
KYRA0100

KYRA Classic™ 500 Non-Fin Stirrups

FOR PATIENTS UP TO 500 LBS. (227 KG.)

The Classic stirrup features an improved easy-squeeze handle and straight rod. The patented design allows maximum rod length to accommodate the widest range of patient heights. This non-finned boot design is often preferred by many robotic surgery teams. Learn about the full line, p. 8.




Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Classic™ 500 Non-Fin Stirrups

KYRA7800

0.75" (19mm) thick pads included. Sold in pairs.
Patient Weight Capacity: 500 lbs. (227 kg)
Stirrup Weight: 14.6 lbs. (6.6 kg)
Lithotomy Range: +87° to -35°
Abduction Range: +25° to -9°
Boot Rotation Angle: +/- 25°



- Non-finned boot design preferred by many robotic surgery teams.
- Easy-squeeze handle with reduced squeeze angle and grip force suits smaller hands.
- Position your shortest to tallest patient.
- Easy-clean silicone straps designed for improved infection control.
- Larger boot with ¾" (19 mm) of pressure-reducing padding for modern patient anatomy.

WHAT'S INCLUDED



Classic 500 Non-Fin Stirrups
w/Secure-Lok™ silicone straps
KYRA7800

WHAT'S REQUIRED



Secure Release™ Clamp*, p. 40
*2 required
KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

WHAT'S OPTIONAL



Purple™ Stirrup Cart, p. 37
Safely stores and transports stirrups.
KYRA1000



Clean™ Catcher Drainage System, p. 34
All-in-one, disposable urology drainage system.
KYRA8452 (Non-sterile)
KYRA8451 (Sterile)



Urology Bands & Drainage Bags, p. 36
KYRA8400 (Band)
KYRA8410 (Non-sterile)
KYRA8411 (Sterile)

POSITION GUIDE

For larger weight capacity, see the
KYRA Classic 800 stirrups →

KYRA COMPLETE Premium Servicing & Warranty

Classic 500 Non-Fin Replacement Pads
Extend the lifetime of your stirrup hardware with replacement stirrup pads.
KYRA6370-P (Secure-Lok Silicone)

KYRAComplete™ Extended Warranty
Cover your hardware and prevent downtime with KYRAComplete premium warranty and service packages. KYRA Stirrups are designed to support common repairs on-site with easy repair kits. Details, p. 80.
KYRA0100

KYRA Classic™ 800 Plus Stirrups

FOR PATIENTS UP TO 800 LBS. (363 KG.)

The Classic stirrup features an improved easy-squeeze handle and straight rod. The patented design allows maximum rod length to accommodate the widest range of patient heights. Extra padding and larger boot size tailored to bariatric patients. Learn about the full line, p. 8.

Fabric-free silicone straps
help improve infection control.



Easy-squeeze Handle
Reduced squeeze angle and grip force suits smaller hands.

Unique Design
Patented design supports safe positioning of the tallest to shortest patients.

Shown with KYRA Clean™ Catcher Drainage System, p. 34

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Classic™ 800 Plus Stirrups

KYRA8850

1" (25mm) thick pads included. Sold in pairs.
Patient Weight Capacity: 800 lbs. (363 kg)
Stirrup Weight: 14.6 lbs. (6.6 kg)
Lithotomy Range: +87° to -35°
Abduction Range: +25° to -9°
Boot Rotation Angle: +/- 25°

- Easy-squeeze handle with reduced squeeze angle and grip force suits smaller hands.
- Position your shortest to tallest patient.
- Easy-clean silicone straps designed for improved infection control.
- Larger boot with 1" (25 mm) of pressure-reducing padding for modern patient anatomy.



WHAT'S INCLUDED



Classic 800 Plus Stirrups
w/Secure-Lok™ silicone straps
KYRA8850

WHAT'S REQUIRED



Secure Release™ Clamp*, p. 40
*2 required
KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

WHAT'S OPTIONAL



Purple™ Stirrup Cart, p. 37
Safely stores and transports stirrups.
KYRA1000



Clean™ Catcher Drainage System, p. 34
All-in-one, disposable urology drainage system.
KYRA8452 (Non-sterile)
KYRA8451 (Sterile)



Urology Bands & Drainage Bags, p. 36
KYRA8400 (Band)
KYRA8410 (Non-sterile)
KYRA8411 (Sterile)

POSITION GUIDE

For larger range of motion, see the KYRA Comfort™ Stirrup Line →

KYRA COMPLETE Premium Servicing & Warranty

Classic 800 Plus Replacement Pads
Extend the lifetime of your stirrup hardware with replacement stirrup pads.
KYRA6380-P (Secure-Lok Silicone)

KYRAComplete™ Extended Warranty
Cover your hardware and prevent downtime with KYRAComplete premium warranty and service packages. KYRA Stirrups are designed to support common repairs on-site with easy repair kits. Details, p. 80.
KYRA0100

KYRA® Comfort™ Stirrups

Maximum range of motion for precision positioning.

Unique curved rod design engineered to support the widest range of lithotomy and abduction—ideal for optimizing positioning for low lithotomy as well as robotics cases.



KYRA COMPLETE Premium Servicing & Warranty

Replacement Pads

Refresh your stirrups without replacing the entire system. When stirrups remain sound, new boot pads restore comfort and protection.


KYRAComplete™ Extended Warranty

Cover your hardware and prevent downtime with KYRAComplete premium warranty and service packages. KYRA Stirrups are designed to support common repairs on-site with easy repair kits. Details, p. 80.

KYRA0100

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

The KYRA Comfort™ Difference




✓ Curved rod provides larger range of motion optimized for low lithotomy and robotics



✓ Free-floating boot design safely maintains proper angle for patient's leg



✓ Flexible lateral fin, increased padding and larger boot for peroneal nerve protection



✓ Easy-clean, fabric-free silicone SecureLok™ straps* help improve infection control and snap in securely with an audible click



✓ Additional 5° of boot rotation internally and externally for proper alignment on wider tables

On-site repairs

✓ Maximize uptime with repair kits, extended warranty and premium servicing

* Certain models available with hook and loop straps.

Find Your Fit

Model Comparison Guide







	Comfort 350 Stirrups	Comfort 500 Stirrups	Comfort 800 Stirrups
PATIENT WEIGHT CAPACITY	350 lbs. (159 kg)	500 lbs. (227 kg)	800 lbs. (363 kg)
SUGGESTED PATIENT BMI	Under 29	29 or greater	29 or greater
STIRRUP WEIGHT	13.2 lbs. (6 kg)	14.6 lbs. (6.6 kg)	14.6 lbs. (6.6 kg)
BOOT PADS THICKNESS	¾" (19 mm) pads	1" (25 mm) pads	1" (25 mm) pads
STRAP OPTIONS	Silicone or Hook & Loop	Silicone or Hook & Loop	Silicone only
LITHOTOMY RANGE	+85° to -55° (22° more than competitive stirrups)		
USE CASE	General lithotomy, low lithotomy and robotics cases.		


Looking for a **fabric-free** stirrup?
Learn about **KYRA Clean™ Stirrups**, p. 26

Complete ordering details →

KYRA Comfort™ 350 Stirrups

FOR PATIENTS UP TO 350 LBS. (159 KG.)

Engineered to support the widest range of lithotomy and abduction—ideal for optimizing positioning, surgical site access and robotics. Learn about the full line, p. 18.



Looking for a stirrup with a squeeze-handle?


Learn about **KYRA Classic™ Stirrups**, p. 8

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Comfort™ 350 Stirrups

KYRA3100

0.75" (19mm) thick pads included. Sold in pairs.
Patient Weight Capacity: 350 lbs. (159 kg)
Stirrup Weight: 13.2 lbs. (5.9 kg)
Lithotomy Range: +85° to -55°
Abduction Range: +25° to -9°
Boot Rotation Angle: +/- 25°



- ✓ Curved rod facilitates maximum range of motion for robotics access, low lithotomy and optimal positioning.
- ✓ Easy-clean, durable silicone straps designed for improved infection control.
- ✓ Larger boot with ¾" (19 mm) of pressure-reducing padding for modern patient anatomy.
- ✓ Accommodates your tallest patients.

WHAT'S INCLUDED



Comfort 350 Stirrups
w/Secure-Lok™ silicone straps
KYRA3100

OR



Comfort 350 Stirrups
w/hook and loop straps
KYRA3000

WHAT'S REQUIRED



Secure Release™ Clamp*, p. 40
*2 required

KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

WHAT'S OPTIONAL



Purple™ Stirrup Cart, p. 37
Safely stores and transports stirrups.
KYRA1000



Clean™ Catcher Drainage System, p. 34
All-in-one, disposable urology drainage system.
KYRA8452 (Non-sterile)
KYRA8451 (Sterile)



Urology Bands & Drainage Bags, p. 36
KYRA8400 (Band)
KYRA8410 (Non-sterile)
KYRA8411 (Sterile)

POSITION GUIDE

For larger weight capacities, see the KYRA Comfort 500 & 800 stirrups →

KYRA COMPLETE Premium Servicing & Warranty

Comfort 350 Replacement Pads
Extend the lifetime of your stirrup hardware with replacement stirrup pads.
KYRA6260-P (Secure-Lok Silicone)
KYRA6210-P (Hook & Loop)

KYRAComplete™ Extended Warranty
Cover your hardware and prevent downtime with KYRAComplete premium warranty and service packages. KYRA Stirrups are designed to support common repairs on-site with easy repair kits. Details, p. 80.
KYRA0100

KYRA Comfort™ 500 Stirrups

FOR PATIENTS UP TO 500 LBS. (227 KG.)

Engineered to support the widest range of lithotomy and abduction—ideal for optimizing positioning, surgical site access and robotics. Learn about the full line, p. 18.



Fabric-free silicone straps
help improve infection control.

Larger Boots
with pressure-reducing padding help protect from nerve damage.

Curved Rod
provides maximum space during low lithotomy and robotics cases.

Low-Torque Twist Handle
reduces strain and fatigue for clinician.


Looking for a stirrup with a squeeze-handle?
Learn about **KYRA Classic™ Stirrups**, p. 8

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Comfort™ 500 Stirrups

KYRA5100

1" (25mm) thick pads included. Sold in pairs.
Patient Weight Capacity: 500 lbs. (227 kg)
Stirrup Weight: 14.6 lbs. (6.6 kg)
Lithotomy Range: +85° to -55°
Abduction Range: +25° to -9°
Boot Rotation Angle: +/- 25°



- Curved rod facilitates maximum range of motion for robotics access, low lithotomy and optimal positioning.
- Easy-clean, durable, fabric-free silicone straps designed for improved infection control.
- Larger boot with 1" (25 mm) of pressure-reducing padding for modern patient anatomy.
- Accommodates your tallest patients.

WHAT'S INCLUDED



Comfort 500 Stirrups
w/Secure-Lok™ silicone straps
KYRA5100

OR



Comfort 500 Stirrups
w/hook and loop straps
KYRA5000

WHAT'S REQUIRED

Secure Release™ Clamp*, p. 40
*2 required
KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

WHAT'S OPTIONAL



Purple™ Stirrup Cart, p. 37
Safely stores and transports stirrups.
KYRA1000



Clean™ Catcher Drainage System, p. 34
All-in-one, disposable urology drainage system.
KYRA8452 (Non-sterile)
KYRA8451 (Sterile)



Urology Bands & Drainage Bags, p. 36
KYRA8400 (Band)
KYRA8410 (Non-sterile)
KYRA8411 (Sterile)

POSITION GUIDE

For larger weight capacity, see the KYRA Comfort 800 stirrups →

KYRA COMPLETE Premium Servicing & Warranty


Comfort 500 Replacement Pads
Extend the lifetime of your stirrup hardware with replacement stirrup pads.
KYRA6250-P (Secure-Lok Silicone)
KYRA6200-P (Hook & Loop)

KYRAComplete™ Extended Warranty
Cover your hardware and prevent downtime with KYRAComplete premium warranty and service packages. KYRA Stirrups are designed to support common repairs on-site with easy repair kits. Details, p. 80.
KYRA0100

KYRA Comfort™ 800 Stirrups

FOR PATIENTS UP TO 800 LBS. (363 KG.)

Extra padding and larger boot size tailored to bariatric patients. Engineered to support the widest range of lithotomy and abduction—ideal for optimizing positioning, surgical site access and robotics. Learn about the full line, p. 18.




Looking for a stirrup with a squeeze-handle?

Learn about **KYRA Classic™ Stirrups**, p. 8

Comfort™ 800 Stirrups

KYRA8100

1" (25mm) thick pads included. Sold in pairs.
Patient Weight Capacity: 800 lbs. (363 kg)
Stirrup Weight: 14.6 lbs. (6.6 kg)
Lithotomy Range: +85° to -55°
Abduction Range: +25° to -9°
Boot Rotation Angle: +/- 25°



- Curved rod facilitates maximum range of motion for robotics access, low lithotomy optimal positioning.
- Easy-clean, durable, fabric-free silicone straps designed for improved infection control.
- Larger boot with 1" (25 mm) of pressure-reducing padding for modern patient anatomy.
- Accommodates your heaviest and tallest patients.

WHAT'S INCLUDED



Comfort 800 Stirrups
w/Secure-Lok™ silicone straps
KYRA8100

WHAT'S REQUIRED



Secure Release™ Clamp*, p. 40
*2 required
KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

WHAT'S OPTIONAL



Purple™ Stirrup Cart, p. 37
Safely stores and transports stirrups.
KYRA1000



Clean™ Catcher Drainage System, p. 34
All-in-one, disposable urology drainage system.
KYRA8452 (Non-sterile)
KYRA8451 (Sterile)



Urology Bands & Drainage Bags, p. 36
KYRA8400 (Band)
KYRA8410 (Non-sterile)
KYRA8411 (Sterile)

POSITION GUIDE

For **fabric-free** stirrups, see the KYRA Clean™ stirrups →

KYRA COMPLETE Premium Servicing & Warranty

Comfort 800 Replacement Pads
Extend the lifetime of your stirrup hardware with replacement stirrup pads.
KYRA6320-P (Secure-Lok Silicone)

KYRAComplete™ Extended Warranty
Cover your hardware and prevent downtime with KYRAComplete premium warranty and service packages. KYRA Stirrups are designed to support common repairs on-site with easy repair kits. Details, p. 80.
KYRA0100

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

KYRA® Clean™ Stirrups

Designed for maximum infection control.

A fully fabric-free stirrup designed for maximum infection control, the Clean Stirrup range delivers all the positioning performance of Comfort Stirrups (p. 18) with easy-clean, coated foam pads that are completely wipe-clean, helping minimize contamination risk.

Easy Disinfection
Seamless, waterproof surfaces wipe clean without crevices or seams.

Fabric-free design
Helps improve infection control with durable, easy-to-clean materials.

Flexible comfort
Soft, coated pads cushion and conform around extremities for secure positioning.



KYRA COMPLETE Premium Servicing & Warranty

Replacement Pads
Refresh your stirrups without replacing the entire system. When stirrups remain sound, new boot pads restore comfort and protection.

KYRAComplete™ Extended Warranty
Cover your hardware and prevent downtime with KYRAComplete premium warranty and service packages. KYRA Stirrups are designed to support common repairs on-site with easy repair kits. Details, p. 80.

KYRA0100

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

The KYRA Clean™ Difference



✓ Antimicrobial, waterproof, coated pads wipe clean quickly between cases



✓ Soft, flexible pads conform around extremities for reliable, secure positioning




✓ Easy-clean, fabric-free silicone Secure-Lok™ straps help improve infection control

On-site repairs

✓ Maximize uptime with repair kits, extended warranty and premium servicing

Same positioning performance as the Comfort™ line, full details, p. 18



✓ Curved rod provides larger range of motion optimized for robotics access
✓ Free-floating boot design safely maintains proper angle for patient's leg
✓ Flexible lateral fin and larger boot for peroneal nerve protection





Find Your Fit

Model Comparison Guide


	Clean 350 Stirrups	Clean 500 Stirrups
PATIENT WEIGHT CAPACITY	350 lbs. (159 kg)	500 lbs. (227 kg)
SUGGESTED PATIENT BMI	Under 29	29 or greater
STIRRUP WEIGHT	13.2 lbs. (5.9 kg)	14.4 lbs. (6.5 kg)
BOOT PADS THICKNESS	¾" (19 mm) pads	1" (25 mm) pads
STRAP OPTIONS	Silicone Secure-Lok™	Silicone Secure-Lok™
LITHOTOMY RANGE	+85° to -55° (22° more than competitive stirrups)	
USE CASE	General lithotomy, low lithotomy and robotics cases designed for maximum infection control.	

Complete ordering details →

KYRA Clean™ 350 Stirrups

FOR PATIENTS UP TO 350 LBS. (159 KG.)

The fabric-free stirrup-the Clean Stirrup range delivers all the positioning performance of Comfort Stirrups (p. 18) with easy-clean, coated foam pads that are completely wipe-clean, helping minimize contamination risk and improve infection control. Learn about the full line, p. 26.



Looking for a stirrup with a squeeze-handle?

Learn about **KYRA Classic™ Stirrups**, p. 8

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Clean™ 350 Stirrups

KYRA3200

0.75" thick (19 mm) antimicrobial-coated foam pads included.
Sold in pairs.
Patient Weight Capacity: 350 lbs. (159 kg)
Stirrup Weight: 13.2 lbs. (5.9 kg)
Lithotomy Range: +85° to -55°
Abduction Range: +25° to -9°
Boot Rotation Angle: +/- 25°



- Fabric-free stirrup with easy-clean, durable pads and silicone straps designed for improved infection control.
- Antimicrobial, vinyl-coated foam pads are removable and fully wipe-clean, minimizing infection risk between cases.
- Curved rod facilitates maximum range of motion for robotics access, low lithotomy optimal positioning.

WHAT'S INCLUDED



Clean™ 350 Stirrups
w/Secure-Lok™ silicone straps
KYRA3200

WHAT'S REQUIRED



Secure Release™ Clamp*, p. 40
*2 required
KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

WHAT'S OPTIONAL



Purple™ Stirrup Cart, p. 37
Safely stores and transports stirrups.
KYRA1000



Clean™ Catcher Drainage System, p. 34
All-in-one, disposable urology drainage system.
KYRA8452 (Non-sterile)
KYRA8451 (Sterile)



Urology Bands & Drainage Bags, p. 36
KYRA8400 (Band)
KYRA8410 (Non-sterile)
KYRA8411 (Sterile)

POSITION GUIDE

For higher weight capacity **fabric-free** stirrups, see the KYRA Clean™ 500 stirrups →

KYRA COMPLETE Premium Servicing & Warranty

Clean™ 350 Replacement Pads
Extend the lifetime of your stirrup hardware with replacement stirrup pads.
KYRA6290-P (Secure-Lok Silicone)

KYRAComplete™ Extended Warranty
Cover your hardware and prevent downtime with KYRAComplete premium warranty and service packages. KYRA Stirrups are designed to support common repairs on-site with easy repair kits. Details, p. 80.
KYRA0100


KYRA Clean™ 500 Stirrups

FOR PATIENTS UP TO 500 LBS. (227 KG.)

The fabric-free stirrup—the Clean Stirrup range delivers all the positioning performance of Comfort Stirrups (p. 18) with easy-clean, coated foam pads that are completely wipe-clean, helping minimize contamination risk and improve infection control. Learn about the full line, p. 26.

Flexible comfort
Soft, coated pads cushion and conform around extremities for secure positioning.

Fabric-free design
Helps improve infection control with durable, easy-clean materials.



Easy Disinfection
Seamless, waterproof surfaces wipe clean without crevices or seams.

Curved Rod
provides maximum space during low lithotomy and robotics cases.


Looking for a stirrup with a squeeze-handle?
Learn about **KYRA Classic™ Stirrups**, p. 8

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Clean™ 500 Stirrups

KYRA5200

1" (25 mm) thick antimicrobial-coated foam pads included.
Sold in pairs.
Patient Weight Capacity: 500 lbs. (227 kg)
Stirrup Weight: 14.4 lbs. (6.5 kg)
Lithotomy Range: +85° to -55°
Abduction Range: +25° to -9°
Boot Rotation Angle: +/- 25°



- Fabric-free stirrup with easy-clean, durable pads and silicone straps designed for improved infection control.
- Antimicrobial, vinyl-coated foam pads are removable and fully wipe-clean, minimizing infection risk between cases.
- Curved rod facilitates maximum range of motion for robotics access, low lithotomy optimal positioning.

WHAT'S INCLUDED



Clean™ 500 Stirrups
w/Secure-Lok™ silicone straps
KYRA5200

WHAT'S REQUIRED



Secure Release™ Clamp*, p. 40
*2 required
KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

WHAT'S OPTIONAL



Purple™ Stirrup Cart, p. 37
Safely stores and transports stirrups.
KYRA1000



Clean™ Catcher Drainage System, p. 34
All-in-one, disposable urology drainage system.
KYRA8452 (Non-sterile)
KYRA8451 (Sterile)



Urology Bands & Drain Bags, p. 36
KYRA8400 (Band)
KYRA8410 (Non-sterile)
KYRA8411 (Sterile)

POSITION GUIDE

For traditional padded stirrups, see the KYRA Comfort™ 500 stirrups p. 23

KYRA COMPLETE Premium Servicing & Warranty


Clean™ 500 Replacement Pads
Extend the lifetime of your stirrup hardware with replacement stirrup pads.
KYRA6270-P (Secure-Lok Silicone)

KYRAComplete™ Extended Warranty
Cover your hardware and prevent downtime with KYRAComplete premium warranty and service packages. KYRA Stirrups are designed to support common repairs on-site with easy repair kits. Details, p. 80.
KYRA0100

KYRA 400 Stirrups

KYRA7100

Pads included • Sold in pairs
Patient Weight Capacity: 400 lbs. (181 kg)
Stirrup Weight: 12.2 lbs. (5.5 kg)
Lithotomy Range: +90˚ to -55˚
Abduction Range: +25˚ to -9˚
Boot Rotation Angle: +/- 25˚



- Non-finned boot design preferred by many robotic surgery teams.
- Easy-clean, durable silicone straps designed for improved infection control.
- ¾" (19 mm) of supportive padding for enhanced patient comfort and support.

POSITION GUIDE

WHAT'S INCLUDED



KYRA 400 Stirrups
w/Secure-Lok™ silicone straps
KYRA7100

WHAT'S REQUIRED



Secure Release™ Clamp*, p. 40
*2 required
KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

WHAT'S OPTIONAL



Purple™ Stirrup Cart, p. 37
Safely stores and transports stirrups.
KYRA1000



KYRA 250 Stirrups

KYRA7810

Pads included • Sold in pairs
Patient Weight Capacity: 250 lbs. (113kg)
Stirrup Weight: 10.2 lbs. (4.6 kg)
Lithotomy Range: +85˚ to -33˚
Abduction Range: +22˚ to -9˚
Boot Rotation Angle: +/- 25˚

- Non-finned boot design sized for pediatric patients aged 6+ and smaller adults up to 250 lbs.
- Easy-clean, durable silicone straps designed for improved infection control.
- ¾" (19 mm) of supportive padding for enhanced patient comfort and support.

POSITION GUIDE

WHAT'S INCLUDED



KYRA 250 Stirrups
w/Secure-Lok™ silicone straps
KYRA7810

WHAT'S REQUIRED



Secure Release™ Clamp*, p. 40
*2 required
KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

WHAT'S OPTIONAL



Pediatric Stirrup Cart, p. 37
Safely stores and transports stirrups.
KYRA1010

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.



Knee Crutch, Lift Assist

KYRA7305

Pads included • Sold in pairs
Patient Weight Capacity: 500 lbs. (227 kg)
Device Weight: 9.6 lbs. (4.3 kg)
Lithotomy Range: +85° to -10°
Abduction Range: +22° to -9°

- Innovative gas spring uses lift-assist to ease lithotomy positioning.
- Improves precision and control with smoother, more gradual adjustments.
- Enhances staff safety by reducing manual strain.
- Included leg strap helps secure patient's legs.

POSITION GUIDE

WHAT'S INCLUDED



Knee Crutch Lift Assist
KYRA7305

WHAT'S REQUIRED



Secure Release™ Clamp*, p. 40
*2 required
KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

WHAT'S OPTIONAL



Pediatric Stirrup Cart, p. 37
Safely stores and transports stirrups.
KYRA1010

Find Your Fit

Model Comparison Guide





	KYRA 250 Stirrups	KYRA 400 Stirrups	Lift-Assist Knee Crutch
PATIENT WEIGHT CAPACITY	250 lbs. (113 kg)	400 lbs. (181 kg)	500 lbs. (227 kg)
DEVICE WEIGHT	10.2 lbs. (4.6 kg)	12.2 lbs. (5.5 kg)	9.6 lbs. (4.6 kg)
BOOT PAD THICKNESS	¾" (19 mm) pads	¾" (19 mm) pads	¾" (19 mm) padding
STRAP OPTIONS	Silicone Secure-Lok™	Silicone Secure-Lok™	Padded safety strap included
SUPPORT TYPE	Flexible boot with non-finned design		Goepel-style
LIFT ASSIST TECHNOLOGY	No	Yes	Yes
LITHOTOMY RANGE	+90˚ to -35˚	+90˚ to -55˚	+85˚ to -10˚
USE CASE	Designed for patients ages 6+, weighing up to 250 lbs.	Non-finned boot for robotics access.	Ideal for urology procedures/ fitment on urology chairs.

For more lithotomy positioning options, see our full line of
Classic, Comfort and Clean stirrups, p. 7-31

KYRA Clean™ Catcher Drainage System

Fully disposable drainage system.

The KYRA® Clean™ Catcher System provides a fully disposable drainage system improving hygiene, reducing fluid exposure, and streamlining setup and disposal.



Allows for proper hip alignment with stirrup when used in conjunction with Clean Catcher mounts.

No more messy clean-ups!

Permanent stirrup mounts saves rail space and ensures optimal patient positioning.

Fully disposable system for improved infection control

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Clean™ Catcher Drainage System

KYRA8452 - Non-Sterile, 10/cs
KYRA8451 - Sterile, 10/cs

Max. Volume Capacity: 2 gallons (8L)



- Fully disposable drainage system designed for maximum infection control.
- No more messy clean-ups or lost urobands.
- Versatile mounting options—attaches to the distal table rail or directly to the stirrup via unique mount system.
- Stable, no-sag performance prevents band droops.

POSITION GUIDE

WHAT'S INCLUDED



Non-sterile Clean™ Catcher Drainage System

KYRA8452

OR



Sterile Clean™ Catcher Drainage System

KYRA8451

WHAT'S OPTIONAL



Clean™ Catcher Stirrup Mounts

Permanent stirrup mounts, details below.

KYRA8460 (Round)
KYRA8461 (Flat)
KYRA8462 (Allen)

Clean™ Catcher Stirrup Mounts

KYRA8460 - Round
KYRA8461 - Flat
KYRA8462 - Allen*

Sold in pairs





- Permanent, dedicated connection point for attachment of the KYRA Clean™ Catcher Drainage System.
- Space-saving adapter allows optimal alignment of stirrup to patient's hip, reducing risk of hip dislocations.
- Factory-installed or field-retrofittable on any KYRA or Allen* stirrups for immediate integration.



Connects directly to stirrup via mount—saving rail space and allowing for optimal patient positioning.

*Allen® is a registered trademark of Baxter.



Shown with Urology Attachment Band, sold below.

Urology Drainage Bags

KYRA8410 - Non-Sterile
KYRA8411 - Sterile

Max. Volume Capacity: 2 gallons (8L)

- Large-capacity drainage bag with integrated buttocks flap for secure, controlled fluid capture.
- Disposable for efficient cleanup and reduced contamination risk.
- Compatible with KYRA Urology Attachment Band (KYRA8400).

POSITION GUIDE

WHAT'S INCLUDED



Non-sterile Urology Drainage Bags
10/case
KYRA8410

OR



Sterile Urology Drainage Bags
10/case
KYRA8411

WHAT'S REQUIRED



Urology Attachment Band, Universal
KYRA8400



Integrated rail clamp universally fits all rail types

Urology Attachment Band, Universal

KYRA8400

Band: 22.5" Diam. x 1.25"H x 0.125" thick (57 cm x 3.2 cm x 0.3 cm)
Integrated Clamp: 1"W x 3.6"H x 1" thick (2.5 cm x 9.1 cm x 2.5 cm)

- Durable, reusable support band engineered to provide stable, reliable support.
- Integrated rail clamps ensure quick, secure attachment to any OR table side rail.
- Compatible with KYRA Urology Drainage bags.

POSITION GUIDE

WHAT'S INCLUDED



Urology Attachment Band, Universal
One included
KYRA8400

WHAT'S REQUIRED



Non-sterile Urology Drainage Bags
10/case
KYRA8410

OR



Sterile Urology Drainage Bags
10/case
KYRA8411

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Purple Stirrup Cart

KYRA1000 - Purple Stirrup Cart
KYRA1010 - Pediatric Stirrup Cart

Cart Height: 40.75" (103.5 cm)
Base Width: 21.625" x 21.625" (54.93 cm x 54.93 cm)
Cart Weight: 35.5 lbs. (16.1 kg)

- Secures stirrups during transport, reducing risk of damage.
- Ergonomic handles for easy maneuverability.
- Integrated clamp storage for added convenience and organization.
- Locking casters prevent unwanted movement.



Clamp storage for accessibility



WHAT'S INCLUDED



Purple Stirrup Cart
KYRA1000

OR



Pediatric Stirrup Cart
KYRA1010

WHAT'S OPTIONAL



Secure Release™ Clamp, p. 40
KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

WHAT'S COMPATIBLE

	Classic™ Stirrups		Comfort™ Stirrups			Clean™ Stirrups		Specialty Stirrups & Leg Holders				
	350	500	350	500	800	350	500	KYRA 400	KYRA 250	Lift Assist Knee Crutch	Simple Stirrup	Well Leg Holder
PURPLE CART	✓	✓	✓	✓	✓	✓	✓	✓	✓	x	x	x
PEDIATRIC CART	x	x	x	x	x	x	x	x	✓	✓	✓	✓



Knee Crutch, Goepel
KYRA7300

Pads included • Sold in pairs
Patient Weight Capacity: 500 lbs. (227 kg)
Device Weight: 6.6 lbs. (3.0 kg)

- Classic Goepel-style knee crutch.
- Included leg strap helps secure patient's legs.
- Removable, easy-to-clean pads help reduce pressure points.

POSITION GUIDE

WHAT'S INCLUDED



Knee Crutch, Goepel

KYRA7300

WHAT'S REQUIRED



SocketLok™ Clamp*, p. 41
*2 required

KYRA4850 (US) KYRA4860 (EU)
KYRA4870 (UK) KYRA4890 (JP)

WHAT'S OPTIONAL



Urology Bands & Drainage Bags, p. 36

KYRA8400K (Band)
KYRA8410 (Non-sterile)
KYRA8411 (Sterile)

Well Leg Holder
KYRA9205

Pad included. Sold individually.
Patient Weight Capacity: 600 lbs. (272 kg)
Device Weight: 6.9 lbs. (3.1 kg)

- Full boot rotation supports and securely positions the non-operative leg out of the surgical working area.
- Easy setup and articulation for precise, efficient positioning.
- Universal design fits left or right leg.

POSITION GUIDE



Full Boot Rotation
Securely position the non-operative leg

WHAT'S INCLUDED



Well Leg Holder
w/Secure-Lok™ silicone straps

KYRA9205

WHAT'S REQUIRED



SocketLok™ Clamp*, p. 41
*One required

KYRA4850 (US) KYRA4860 (EU)
KYRA4870 (UK) KYRA4890 (JP)

WHAT'S OPTIONAL



Pediatric Stirrup Cart, p. 37

Safely stores and transports stirrups.


KYRA1010

Simple Stirrups
KYRA7860

Pads included • Sold in pairs.
Patient Weight Capacity: 600 lbs. (272 kg)
Stirrup Weight: 7 lbs. (3.18 kg)

- Non-finned boot design with simple working mechanism.
- Easy-clean, durable silicone straps designed for improved infection control.
- ¾" (19 mm) of supportive padding for enhanced patient comfort and support.

POSITION GUIDE



WHAT'S INCLUDED



Simple Stirrups
w/Secure-Lok™ silicone straps

KYRA7860

WHAT'S REQUIRED



SocketLok™ Clamp*, p. 41
*2 required

KYRA4850 (US) KYRA4860 (EU)
KYRA4870 (UK) KYRA4890 (JP)

WHAT'S OPTIONAL



Pediatric Stirrup Cart, p. 37

Safely stores and transports stirrups.

KYRA1010

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Flat Blade Clamps—1" (2.5cm) x 0.25" (0.6 cm)



Secure Release™ Clamp
KYRA4100

Sold individually.
Accepts flat blades, 1" x 0.25" (2.5 x 0.6 cm).
Patient Weight Capacity: 800 lbs. (363 kg)
Device Weight: 0.86 lbs. (0.39 kg)

- ✓ Secure Release™ button prevents false rail engagement & simplifies removal.
- ✓ Mounts and locks anywhere on a standard OR table rail for maximum convenience.
- ✓ Clearly labeled for rail compatibility (US, UK, EU, JP, DY).

Specify Rail Type

US	KYRA4100	UK	KYRA4300	JP	KYRA4500
EU	KYRA4200	DY	KYRA4400		



BladeLok™ Clamp
KYRA4650

Sold individually.
Accepts flat blades, 1" x 0.25" (2.5 x 0.6 cm).
Patient Weight Capacity: 800 lbs. (363 kg)
Device Weight: 0.65 lbs. (0.29 kg)

- ✓ Quick-release clamp mechanism safely mounts and locks anywhere on the table rail for convenience.
- ✓ Clamp locks independently to rail allowing for secure positioning.

Specify Rail Type

US	KYRA4650	UK	KYRA4670
EU	KYRA4660	JP	KYRA4690



Simple Blade Clamp
KYRA4700

Sold individually.
Accepts flat blades, 1" x 0.25" (2.5 x 0.6 cm).
Patient Weight Capacity: 800 lbs. (363 kg)
Device Weight: 0.43 lbs. (0.20 kg)

- ✓ Secures flat post accessories at table rails ends or rail notches.
- ✓ Simple screw system for quick, reliable setup.
- ✓ Cost-effective option.

Specify Rail Type

US	KYRA4700	EU/UK	KYRA4710
DY	KYRA4730	JP	KYRA4740

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Round Post Clamps & Side Rail Extenders



SocketLok™ Clamp
KYRA4850

Sold individually.
Accepts round posts up to 0.6875" (1.8 cm).
Patient Weight Capacity: 800 lbs. (363 kg)
Device Weight: 1.26 lbs. (0.57 kg)

- ✓ Quick-release clamp mechanism safely mounts and locks anywhere on the table rail for convenience.
- ✓ Clamp locks independently to rail allowing for secure positioning of round post accessories.

Specify Rail Type

US	KYRA4850	UK	KYRA4870
EU	KYRA4860	JP	KYRA4890



Side Rail Extender
KYRA2305

Clamp included. Sold individually.
Patient Weight Capacity: 800 lbs. (363 kg)
Dimensions: 8" x 3" x 3" (20 cm x 7.6 cm x 7.6 cm)

- ✓ Extend any surgical table rail by 6" (15.2 cm) for added accessory flexibility.
- ✓ Secure, and simple to attach—no tools required.


Specify Rail Type

US	KYRA2305	EU	KYRA2306

Table Compatibility Guide

Standard O.R. Rail Dimensions

Rail Type	Height		Width	
AMERICAN (US)	1.125" (1⅛")	28.6 mm	0.375" (⅜")	9.5 mm
ESCHMANN (UK)	1.25" (1¼")	31.75 mm	0.25 (¼")	6.35 mm
EUROPEAN (EU)	0.98"	25.0 mm	0.40" (~⅜")	10.0 mm
JAPANESE (JP)	1.25" (1¼")	32.0 mm	0.35" (⅜")	9.0 mm
DENYERS (DY)	1.5" (1½")	38.1 mm	0.25 (¼")	6.35 mm



Alternate order numbers are provided for country-specific rail versions.



Flexible Anesthesia Screen
KYRA2501

Sold individually.
Dimensions: 66.5"L x 1"W (169 cm x 2.5 cm) • Screen Loop: 60"L (152 cm)
Device Weight: 2 lbs. (0.9 kg)

- Coated gooseneck band provides flexible positioning for anesthesia access.
- Helps maintain a consistent separation between the anesthesia and surgical fields.
- Durable, lightweight construction supports efficient setup and easy cleaning.

POSITION GUIDE

WHAT'S INCLUDED



Flexible Anesthesia Screen
KYRA2501

WHAT'S REQUIRED



SocketLok™ Clamp*, p. 41
*2 required
KYRA4850 (US) KYRA4860 (EU)
KYRA4870 (UK) KYRA4890 (JP)

Table Width Extender
KYRA2320

Pad included • Sold individually • US Rail only
Patient Weight Capacity: 800 lbs. (363 kg)
Dimensions: 20"L x 4"W (50.8 cm x 10 cm)



- Extends OR table width by 4" (10.2 cm) to support larger patient surface area.
- Integrated side rail provides secure accessory attachment.
- Includes matching pad for added patient comfort and support.

POSITION GUIDE

WHAT'S INCLUDED




Table Width Extender
KYRA2320

+



Pad for Table Width Extender, 2" (5 cm) thick
KYRA2321

WHAT'S OPTIONAL



Arm Protection Shield, p. 43
w/padded cover
KYRA2330

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Arm Protection Shield
KYRA2330

Pad included • Sold individually.
Dimensions: 11" x 16" x 0.9" (28 cm x 40 cm x 23 cm)



- Provides stable containment to protect the ulnar nerve and maintain arm position.
- Black padded cover slips over the clear shield, adding protection while preserving access for anesthesia and surgical teams.
- Durable, lightweight design allows easy positioning and efficient cleaning.

POSITION GUIDE

WHAT'S INCLUDED



Arm Protection Shield
w/padded cover
KYRA2330

WHAT'S RELATED



Clean™ Catcher Drainage System, p. 34
All-in-one, disposable urology drainage system.
KYRA8452 (Non-sterile)
KYRA8451 (Sterile)



Urology Bands & Drain Bags, p. 36
KYRA8400 (Band)
KYRA8410 (Non-sterile)
KYRA8411 (Sterile)

Bariatric Armboard
KYRA2060

Pad included • Sold individually • US Rail only
Patient Weight Capacity: 800 lbs. (363 kg)
Dimensions: 26"L x 8"W (66 cm x 20 cm)
Armboard Weight: 6 lbs. (2.7 kg)



2.5" (6.4cm) wider than standard armboards

- Wider version of our trigger-style armboard, ideal for bariatric patients.
- Integrated clamp snap locks onto rail at any position.
- Single-handed adjustment simplifies set-up.
- Eliminates pinch hazards and gear grinding.
- Durable, radiolucent construction.
- Includes 2" (5cm) thick pad for patient comfort.

POSITION GUIDE

WHAT'S INCLUDED



Bariatric Armboard
Pad included
KYRA2060

WHAT'S OPTIONAL



Soft Disposable Armboard Straps, p.79
12/pack
KYRA8310

Trigger Armboard
KYRA2000

Pad not included. Sold individually.
Dimensions: 26"L x 5.5"W (66 cm x 14 cm)
Armboard Weight: 4.0 lbs. (1.8 kg)
Patient Weight Capacity: 500 lbs. (227 kg)

- ✓ Integrated clamp locks onto rail at any position.
- ✓ Waterfall design prevents trapped liquids and contamination.
- ✓ No exposed teeth—eliminates pinch hazards and gear grinding
- ✓ Durable, radiolucent construction.



Pivots 180°, shown with Comfort Armboard Pad

Specify Rail Type

US KYRA2000 UK KYRA2006 JP KYRA2008
EU KYRA2005 DY KYRA2007

POSITION GUIDE

WHAT'S INCLUDED



Trigger Armboard
Available in five rail types, see above

WHAT'S REQUIRED



Comfort™ Armboard Pad
Traditional pad in three thicknesses: 1" (2.5 cm), 2" (5 cm) or 3" (7.6 cm)

OR



Clamshell Armboard Pad
Pads with clamshell strap to encapsulate arm, available in 1" (2.5 cm), 2" (5 cm) or 3" (7.6 cm) thick versions

See complete **armboard pad** details, next page →

How it works:



One-handed adjustment

Molded construction helps reduce imaging artifacts

Integrated snap-lock clamp

No metal for radiolucent imaging

For greater flexibility and multi-axis control, see our
Swivel & Articulating Armboards
p. 46



Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

CHOOSE YOUR ARMBOARD PAD



↑ Clamshell strap adds security and protection



Shown on Trigger Armboard, sold separately

Clamshell Armboard Pad
KYRA2050

Sold individually.
Fits armboards: 26"L x 5.5"W (66 cm x 14 cm)
Clamshell Padded Area: 10"W (25.4 cm)
Clamshell Pad Weight:
0.95 lbs. (0.43 kg) for KYRA2050/2051 • 1.0 lbs. (0.45 kg) for KYRA2052

- ✓ Pressure-relieving foam with integrated clamshell strap safely encapsulates arm.
- ✓ Available in three thicknesses—1" (2.5 cm), 2" (5 cm) or 3" (7.6 cm).

Specify Pad Size

1" (2.5 cm) KYRA2050 2" (5 cm) KYRA2051 3" (7.6 cm) KYRA2052

WHAT'S INCLUDED



Clamshell Armboard Pad
1"-3" thick pads with clamshell strap to encapsulate arm.

WHAT'S REQUIRED



Trigger Armboard
Available in five rail types.

See left page.

Comfort™ Armboard Pad
KYRA2009

Sold individually.
Fits armboards: 26"L x 5.5"W (66 cm x 14 cm)
Standard Pad Weight:
0.8 lbs. (0.36 kg) for KYRA2009/2010 • 0.9 lbs. (0.41 kg) for KYRA2015



- ✓ Traditional armboard pads with soft, pressure-relieving foam.
- ✓ Pair of Soft Armboard Straps included
- ✓ Available in three thicknesses—1" (2.5 cm), 2" (5 cm) or 3" (7.6 cm).

Specify Pad Thickness

1" (2.5 cm) KYRA2009 2" (5 cm) KYRA2010 3" (7.6 cm) KYRA2015



Shown with Trigger Armboard, sold separately

WHAT'S INCLUDED



Comfort™ Armboard Pad
Traditional pad in three thicknesses.

WHAT'S REQUIRED



Trigger Armboard
Available in five rail types.

See left page.

WHAT'S OPTIONAL



Soft Disposable Armboard Straps, p. 79
12/pack
KYRA8310

Swivel Armboard

KYRA2304 - Small
KYRA2300 - Large



Same armboard support.
Multiple post reach options.


1" thick Clamshell Pad included. Sold individually.
Overall Dimensions: 16"L x 6"H x 23"W (40.6 cm x 15.2 cm x 58.4 cm)
Small Armboard Post Size: 16" x 6.5" x 5/8" (40.6 cm x 16.5 cm x 1.6 cm)
Large Armboard Post Size: 18.2" x 12.5" x 5/8" (46.2 cm x 31.7 cm x 1.6 cm)
Armboard Weight: Small: 6.0 lbs. (2.7 kg) • Large: 6.1 lbs. (28 kg)
Patient Weight Capacity: 600 lbs. (272 kg)

- ✓ Two swivel post sizes for optimized positioning range.
- ✓ Multi-task armboard, typically used for lateral decubitus and beachchair positions.
- ✓ 1" (25 mm) thick Clamshell Pad included for comfort and protection.

POSITION GUIDE




WHAT'S INCLUDED



Swivel Armboard, Small
One armboard with clamshell pad
KYRA2304

OR



Swivel Armboard, Large
One armboard with clamshell pad
KYRA2300

WHAT'S REQUIRED



SocketLok™ Clamp*, p. 41
*One required
KYRA4850 (US) KYRA4860 (EU)
KYRA4870 (UK) KYRA4890 (JP)

Connect™ Swivel Armboard

KYRA2302

1" thick Clamshell Pad included. Sold individually.
Overall Dimensions: 16" x 7.5" x 4.6" (40.6 cm x 19.05 cm x 11.68 cm)
Armboard Weight: 7.0 lbs. (3.2 kg)
Patient Weight Capacity: 600 lbs. (272 kg)

- ✓ Designed for optimal positioning without C-arm interference.
- ✓ Lightweight, low-profile design.
- ✓ Clamshell Pad included for comfort and protection.
- ✓ Integrated SocketLok™ clamp for quick table attachment.



Low-profile post enables use during proning without interference with imaging

Specify Rail Type

US KYRA2302-US EU KYRA2302-EU
UK KYRA2302-UK JP KYRA2302-JP

POSITION GUIDE

WHAT'S INCLUDED



Connect Swivel Armboard
One armboard with clamshell pad

+



SocketLok™ Clamp
Permanently attached to armboard

WHAT'S OPTIONAL



Connect™ Rail Clamp, 6", p. 72
KYRA9525 (US)
KYRA9527 (EU)
KYRA9529 (JP)

Find Your Fit

Model Comparison Guide





	Swivel, Small	Swivel, Large	Connect Swivel
POSITIONING RANGE	Universal	Universal	Prone
ARMBOARD PAD STYLE	Clamshell Pad, 1" (25mm) thick	Clamshell Pad, 1" (25mm) thick	Clamshell Pad, 1" (25mm) thick
PAD INCLUDED	Yes	Yes	Yes
POST SIZE	16" H x 6.5"W (40.6 cm x 16.5 cm)	18.2"H x 12.5"W (46.2 cm x 31.7 cm)	n/a
INTEGRATED CLAMP	No	No	Yes
USE CASE	Basic support, e.g. well arm in beachchair	Multi-task flexibility, e.g. lateral positioning	Prone positioning for use with C-Arm imaging

Need quick, single-handed positioning?
See our traditional Armboards
p. 44



Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.



Arthroscopic Leg Holder

KYRA9200


Overall Dimensions: 12"L x 2.5"W x 7.5"H (30 cm x 6.4 cm x 19 cm)
Patient Weight Capacity: 500 lbs (227 kg)
Device Weight: 4 lbs. (1.8 kg)

- Designed to stabilize the leg for proper alignment during knee arthroscopy.
- Easy-to-adjust design supports efficient OR setup and positioning.
- Holds leg up to 35" (89 cm) in circumference.



POSITION GUIDE

WHAT'S INCLUDED



Arthroscopic Leg Holder
One included
KYRA9200

WHAT'S REQUIRED



SocketLok™ Clamp*, p. 41
*One required
KYRA4850 (US) KYRA4860 (EU)
KYRA4870 (UK) KYRA4890 (JP)



Disposable Pads
12 per pack, see details below.
KYRA8501

Disposable Pad for Arthroscopic Leg Holder

KYRA8501 - 12 per pack

Dimensions: 37"L x 3.5"W x 1.25"H (94 cm x 8.9 cm x 3.2 cm)

- Non-sterile disposable pad helps to protect from pressure and skin shear.
- Provide added comfort and helps reduce cross-contamination risk.
- For use with the Arthroscopic Leg Holder.



WHAT'S INCLUDED



Disposable Pads
12 per pack
KYRA8501

WHAT'S REQUIRED



Arthroscopic Leg Holder
One included
KYRA9200

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Patient Transfer Boards

Roller Free • Anti-static • Anti-bacterial • Anti-fungal • Latex free • Flame retardant • No metal for MRI use

Not available in CE countries.

- Enables safe, "no-lift" lateral transfers to help reduce caregiver strain.
- Rollboard design requires only a fraction of the patient's body weight to move—typically just 15%.
- Allows a single staff member to transfer average-size patients with minimal effort.
- Smooth, continuous surface helps support patient comfort during transfer.



Cover moves with the patient and rotates around the board

Long/Wide Board

KYRA2401


Board Size: 48"L x 20"W (120 cm x 50 cm)
Board Weight: 3.3 lbs. (1.5 kg)



Short/Wide Board

KYRA2404

Board Size: 40"L x 20"W (100 cm x 50 cm)
Board Weight: 3.3 lbs. (1.5 kg)



Short/Narrow Board

KYRA2403

Board Size: 40"L x 16"W (100 cm x 40 cm)
Board Weight: 2.7 lbs. (1.3 kg)



Nissen Straps

KYRA2603 - Standard
KYRA2602 - Bariatric

Standard: 52"L x 2"W (132 cm x 5 cm) • Pad: 4"W x 16"L x 1.125"H (10.1 cm x 40.6 cm x 2.9 cm) • Weight: 3 lbs. (1.4 kg) • Patient Weight Capacity: 500 lbs. (227 kg)

Bariatric: 73"L x 2"W (185 cm x 5 cm) • 7"W x 32"L x 1.625"H (17.8 cm x 81.3 cm x 4.1 cm) • Weight: 4 lbs. (1.8 kg) • Patient Weight Capacity: 800 lbs. (363 kg)

- Helps prevent patient slipping during Steep Reverse Trendelenburg positioning.
- Easily mounts to any standard OR table rail for quick setup—no carabiners needed.
- Available in two models to support a range of patient sizes.

POSITION GUIDE



Bariatric size shown





Flat blade mount secures strap into side rail clamp

WHAT'S INCLUDED



Standard Nissen Strap
One pair included
KYRA2603

OR



Bariatric Nissen Strap
One pair included
KYRA2602

WHAT'S REQUIRED



Secure Release™ Clamp*, p. 40
*2 required
KYRA4100 (US) KYRA4200 (EU)
KYRA4300 (UK) KYRA4400 (DY)
KYRA4500 (JP)

Secondary Shoulder Supports

KYRA2340

Pads included • Sold in pairs
Overall Dimensions: 9"L x 7"W x 14"H (23 cm x 18 cm x 35 cm)
Patient Weight Capacity: 600 lbs. (272 kg)
Device Weight: 5 lbs. (2.3 kg)

- Provides non-contact shoulder support for extreme Trendelenburg positioning.
- Flexible swivel adjustment enables precise positioning.
- Pads included to help ensure proper pressure distribution.

POSITION GUIDE



WHAT'S INCLUDED



Secondary Shoulder Supports
KYRA2340

WHAT'S REQUIRED



SocketLok™ Clamp*, p. 41
*2 required
KYRA4850 (US) KYRA4860 (EU)
KYRA4870 (UK) KYRA4890 (JP)

WHAT'S OPTIONAL



Arm Protection Shield, p. 43
w/padded cover
KYRA2330



Patient Safety Strap, Medium

KYRA2601

Overall Dimensions: 96"L x 5"W x .5"H (243 cm x 12.7 cm x 1.2 cm)
Device Weight: 2 lbs. (0.9 kg)
Patient Weight Capacity: 500 lbs. (227 kg)
Pad center section: 18"L x 5"W x 0.5" (46 cm x 12.7 cm x 1.3 cm)
Straps are adjustable from 18" to 45" (46 cm to 114 cm)

- Padded safety strap securely holds the patient's torso or legs in place.
- Fast, secure attachment to the surgical table side rail.
- Universal compatibility with all surgical tables and extensions with side rails.

POSITION GUIDE

WHAT'S INCLUDED



Patient Safety Strap, Medium
One included
KYRA2601

WHAT'S RELATED



Trigger Armboard, p. 44
Integrated clamp snap locks onto rail at any position. Five rail types.
KYRA2000



Flexible Anesthesia Screen, p. 42
Coated gooseneck band for flexible positioning
KYRA2501

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.


KYRA Lateral Positioning Set

Secure Support. Safer Lateral Positioning.

The KYRA Lateral Positioning Set delivers stable, streamlined support designed to optimize surgical access while protecting patient safety. Engineered for efficiency and ease of use, KYRA simplifies lateral positioning without compromising control.



How it works:





Lateral Positioning Set
KYRA9100K

3-pc set includes two rectangular and one circular positioner.
Patient Weight Capacity: 600 lbs. (272 kg)
Dimensions: 17"L x 17"H x 4"W (43 cm x 43 cm x 10 cm)
Device Weight: Circular 4.0 lbs. (1.8 kg) • Rectangular 4.1 lbs. (1.86 kg)

- Provides secure, comfortable support for procedures requiring total lateral stability.
- Fully adjustable to any angle with swivel heads that conform to patient anatomy.
- Heads lock securely for precise positioning.
- Silicone covers protect hardware.

POSITION GUIDE

WHAT'S INCLUDED



Rectangular Lateral Positioner
Two included
KYRA9120

+



Circular Lateral Positioner
One included
KYRA9110

WHAT'S REQUIRED



SocketLok™ Clamp*, p. 41
*3 required
KYRA4850 (US) KYRA4860 (EU)
KYRA4870 (UK) KYRA4890 (JP)

WHAT'S OPTIONAL



Swivel Armboard, p. 46
Multi-task armboard with Clamshell pad
KYRA2304 (Small)
KYRA2300 (Large)



Arm Protection Shield, p. 43
w/padded cover
KYRA2330



Flexible Anesthesia Screen, p. 42
Coated gooseneck band for flexible positioning
KYRA2501

Silicone cover protects hardware from fluid and debris.





Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

KYRA Curve™ Spine Frame System

Radiolucent system with multi-platform compatibility.

FOR PATIENTS UP TO 600 LBS. (272 KG.)

The KYRA Curve Spine Frame provides radiolucent, optimal positioning across a wide range of spine procedures. Adjustable height and width enhances intervertebral access and reduces abdominal pressure for a range of patient anatomies. Available in two sizes.



Two Sizes
offer optimal patient positioning for a range of patient anatomies.

How it works:



Medial to lateral adjustment to accommodate patient “belly-hang”.




Durable silicone straps designed for improved infection control.

Curve Cart holds all components securely in one organized, mobile cart.



Locking wheels keep cart stationary



Two sizes offer optimal patient positioning for a range of patient anatomies.

† Compatible with Allen/Jackson-style operating tables with 1.25” (3.2cm) wide and 1.5” (3.8 cm) tall rail sizes.
Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Curve™ Spine Frame System

KYRA9000K - Curve™ Complete Kit
KYRA9005K - Mini Curve™ Complete Kit

Pads included.
Dimensions: Standard - 30”L x 17.5”W x 10”H (76.2 cm x 44.4 cm x 25.4 cm) • Mini - 24”L x 20”W x 7¾”H (61 cm x 51 cm x 19.7 cm)
Patient Weight Capacity: Standard - 600 lbs. (272 kg) • Mini - 500 lbs. (227 kg)
Device Weight: Standard - 24.4 lbs. (11 kg) • Mini - 22 lbs. (10 kg)



Standard Size

Mini Size
(30% shorter)

- Ideal for spine procedures requiring reverse lordosis including laminectomies, decompressions, fusions and discectomies.
- Intraoperative height adjustment for optimal surgical site exposure.
- C-arm compatible and fits both dual-rail spine systems and standard OR tables.
- Easy-clean silicone straps designed for improved infection control.
- Mini version offers a 30% shorter frame to ensure optimal positioning for patients under 5'9".

DUAL RAIL COMPATIBLE† POSITION GUIDE

WHAT'S INCLUDED

Each system includes one Spine Frame (Standard or Mini) with pads, crank handle, strap buckle and storage cart.



Curve Spine Frame
For patients 5'7" and above, up to 600 lbs. (272 kg).
KYRA9000

OR



Mini Curve Spine Frame
For patients 5'9" and under, up to 500 lbs. (227 kg).
KYRA9005



Curve Crank Handle
One rotational handle for controlled, secure lordosis positioning.
KYRA9010

+



Curve Spine Frame Cart
Safely stores and transports spine frame and components.
KYRA1100 (Standard)
KYRA1110 (Mini)

+



Curve Strap Buckle (Pair)
For use with tables without side rails. One pair included.
KYRA9030

WHAT'S OPTIONAL

Curve™ Comfort Covers

KYRA9020S - Standard size
KYRA9024 - Mini size

6 sets per pack

- Layer of supportive foam absorbs excess moisture, dissipates heat and reduces risk of skin shear.
- Single-use for improved infection control.



NEW!

KYRA Connect™ SPINE SYSTEM

A space-saving mobile unit that transforms your OR table into an operative spine & imaging platform



Shown with the Connect™ Prone Support Set and KYRA Prone Head Support System.



Shown with the Curve™ Spine Frame and KYRA Prone Head Support System



Shown with the Connect™ Supine Top

Learn more about the KYRA Connect™ SPINE SYSTEM →

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

KYRA Connect™ SPINE SYSTEM

Bringing unmatched functionality and flexibility to the OR table.

FOR PATIENTS UP TO 500 LBS. (227 KG.)

The KYRA Connect Spine System transforms most standard OR tables into a versatile spine and imaging platform—without the cost or storage footprint of a traditional spine table.

The functionality of a "Jackson-style" table with the storage footprint of an IV pole.



POSITION GUIDE

The KYRA Difference—intuitive, multi-patented design



Table Mounts securely attach Connect rails to existing OR table



Patented 5-axis flex joint supports range of movement of the table without restriction



Integrated bubble level to assist in leveling table



Lightweight, mobile system with footprint of an IV pole easily transports and stores accessories



Auto-locking casters for ease of use and enhanced security



Simple handheld remote controls locking, storage and height adjustments

Powered column provides optimal surgical and imaging height



Supports prone, lateral & supine positioning. Details, p. 70.

Ordering details for the KYRA Connect™ SPINE SYSTEM →

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Connect™ Spine System

KYRA9500

Patient Weight Capacity: 500 lbs. (227 kg)
Device Weight: 100 lbs. (50 kg)
Frame Dimensions (set-up): 69.5"L x 24.4"W x 35"H (176.5 cm x 62.2 cm x 88.9 cm)
Frame Dimensions (parked): 29.1"L x 24.4"W x 62"H (73.9 cm x 62.2 cm x 157.5 cm)

- Transform most standard OR tables into a versatile, cost-effective spine and imaging platform.
- Designed for simple attachment to OR table.
- Storage footprint of an IV pole.
- Compatible with most dual-rail spine accessories.
- Provides full O-arm/C-arm access.
- On-board battery allows for cord-free operation.

POSITION GUIDE

Ideal for overflow cases and increasing capacity



WHAT'S INCLUDED



Connect Spine System Frame

KYRA9500

WHAT'S REQUIRED



Power Cord*
*One required
KYRA9540 (US)
KYRA9541 (EU)
KYRA9542 (JP)
KYRA9543 (UK)



Table Mount Clamps*, p. 73
KYRA9511 (US)
KYRA9512 (EU)
KYRA9515 (JP)

A flexible, configurable system

Customize your KYRA Connect setup by selecting your components, see opposite page.

Prone Supports—ideal for procedures that require flexing or tilting the table



Curve Spine Frame—ideal for laminectomy and discectomy procedures



Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

FIND YOUR FIT

BUILD YOUR KIT

Mix and match accessories to meet your procedural requirements.

Head Supports

Reusable and disposable options for head and neck support.

p. 62



Torso Supports

Prone or supine supports for various surgical requirements.

p. 66



Arm Supports

Multi-axis arm supports allow for optimal positioning.

p. 68



Clamps & Accessories

Expand functionality with clamps and accessories.

p. 72



Disposables

Custom disposables for improved patient care, equipment protection and infection control.

p. 74



connect™ SPINE SYSTEM

Head Supports

KYRA Prone Head Support Systems

Optimal positioning and protection for the prone patient.

FOR PATIENTS UP TO 600 LBS. (272 KG.)

The Prone Head Support System evenly distributes pressure across the patient's face while in the prone position, minimizing risk of pressure injury. The height adjustable design promotes safe, neutral alignment of the patient's head and neck and mirror allows observation of the patient's eyes and tubing.



Patented safety rim



No sharp edges

← Molded design prevents contact with helmet

Shown with Head Plate for use on dual-rail systems, p. 64

Increased height variability compared to competitor models



Easy-turn knobs for simple height adjustment to achieve neutral neckline

† Compatible with Allen/Jackson-style operating tables with 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall rail sizes.
Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Prone Head Support System

KYRA9071K

Dimensions: Helmet: 11.5"L x 9.5"W x 6" H (29.2 cm x 24.1cm x 15.2 cm)
Mirror: 8.5"L x 8.5"W x 2.25"H (21.6 cm x 21.6 cm x 5.7 cm)
Product Weight: Helmet: 1.0 lbs. (0.5 kg) • Mirror: 1.0 lbs. (0.5 kg)

- ✓ Unique molded design with safety rim; prevents contact with helmet, improves protection during transfer.
- ✓ Ergonomic, easy-turn knobs alter height to achieve a neutral neckline.
- ✓ Evenly distributes pressure across the patient's face to reduce the risk of pressure-related injury.
- ✓ Adult and child cushions available to suit a wide range of patients.

DUAL RAIL COMPATIBLE†

POSITION GUIDE



↖ Mirror allows for visualization of eyes and ET tubing

WHAT'S INCLUDED



Prone Helmet
KYRA9071

+



Prone Adjustable Mirror
KYRA9070

WHAT'S REQUIRED



Prone Face Cushion
Ten per pack

KYRA9080 (Adult)
KYRA9083 (Child)

Prone Head Support System, Chinless

KYRA9072K

Dimensions: Helmet: 11.5"L x 9.5"W x 6" H (29.2 cm x 24.1cm x 15.2 cm)
Mirror: 8.5"L x 8.5"W x 2.25"H (21.6 cm x 21.6 cm x 5.7 cm)
Product Weight: Helmet 1.0 lbs. (0.5 kg) • Mirror: 1.0 lbs. (0.5 kg)

Same positioning benefits of the Prone Head Support System plus:

- ✓ Chin-free profile preserves intubation access.
- ✓ Saves time by eliminating the need to remove and reinsert the airway during repositioning.
- ✓ Improves patient safety by reducing contamination risk and avoiding unnecessary airway trauma.

DUAL RAIL COMPATIBLE†

POSITION GUIDE



↖ Chinless design preserves intubation access

Adult size only

WHAT'S INCLUDED



Prone Helmet, Chinless
KYRA9072

+



Prone Adjustable Mirror
KYRA9070

WHAT'S REQUIRED



Prone Face Cushion, Chinless
Ten per pack

KYRA9081 (Adult)

Head Supports

Connect™ Prone Head Plate

KYRA9050

Dimensions: 19”L x 12.6”W x 2.5”H (48.3 cm x 32 cm x 6.4 cm)
Product Weight: 2.3 lbs. (1 kg)

- Easy-release locking clips for easy installation, adjustment and removal.
- Fully radiolucent for clear, unobstructed imaging.
- Optional pivoting mirror allows observation of patients position, eyes and tubing.
- Pivoting mirror features adjustable screws for secure positioning.

DUAL RAIL COMPATIBLE†

POSITION GUIDE



WHAT'S INCLUDED


Connect Prone Head Plate
KYRA9050

+



WHAT'S OPTIONAL

Connect Pivoting Mirror
KYRA9067





Shown with Prone Head Support System for maximum support and stability

ProneAlone™ Head Support

KYRA9082 - 8 per pack

Dimensions: 11.3”L x 9.7”W x 6.5” H (28.7 cm x 24.6 cm x 16.5 cm)
Product Weight: <1 lb. (<0.5kg)

- All-in-one standalone head support solution.
- Molded design for improved support, stability and alignment.
- Fully disposable cushion and mirror for improved infection control.
- Mirror included for facial monitoring.

POSITION GUIDE



WHAT'S INCLUDED



ProneAlone Head Support
8 per pack
KYRA9082

WHAT'S OPTIONAL



Connect Head Plate, p. 64
For use on dual rail systems.
KYRA9050



WHAT'S REQUIRED

HEAD PROTECTION









	Prone Head System	Prone Head System, Chinless	ProneAlone	Foam Head Block
REUSABLE/ DISPOSABLE	Reusable Helmet with Disposable Cushion		Disposable Foam	Disposable Foam
QUANTITY	System includes one Helmet & Mirror		8 per pack	6 per pack
REQUIRED	KYRA9080—Adult Cushion KYRA9083—Child Cushion (10 per pack)	KYRA9081 —Chinless Cushion (10 per pack)	—	—
MIRROR	Mirror included	Mirror included	Disposable mirror	Use with KYRA9050/9067
PATIENT SIZE	Adult and Child Sizes	Adult	Adult	Universal
USE CASE	Maximum stability and support for 8 hours.	Maximum stability and support for 8 hours—chinless design preserves intubation during transfer.	Fully disposable, standalone solution with mirror included—for use up to 6 hours.	Fully disposable cushion for cases up to 2 hours.
FULL DETAILS	See p. 63	See p. 63	See opposite page	See opposite page

† Compatible with Allen/Jackson-style operating tables with 1.25” (3.2cm) wide and 1.5” (3.8 cm) tall rail sizes.
Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.



Foam Head Block

KYRA9022- 6 per pack

Dimensions: 10.5”L x 9.75”W x 6”H (26.7 cm x 24.7 cm x 15.4 cm)
Product Weight: <1 lb. (<0.5kg)

- Supports the head and neck during prone positioning.
- Disposable for improved infection control.
- Slotted design to accommodate ventilation tubing.

POSITION GUIDE

WHAT'S INCLUDED



Foam Head Block
6 per pack
KYRA9022

WHAT'S OPTIONAL



Connect Prone Head Plate System, p. 64
Includes pivoting mirror for patient observation.
KYRA9050K

Torso Supports

Connect™ Prone Support Set

KYRA9052K - Standard
KYRA9053K - Large

Patient Weight Capacity: 600 lbs. (272 kg)
Chest Support: 19"L x 7.5"W x 6.5"H (48.3 cm x 19.1 cm x 16.5 cm) • 1.8 lbs (0.8 kg)
Hip Support: 9"L x 5.75"W x 8.3"H (22.9 cm x 14.6 cm x 21.1 cm) • 1.9 lbs (0.9 kg)
Large Hip Support: 11.25"L x 5.75"W x 4"H (28.6 cm x 14.6 cm x 10.2 cm) • 2.4 lbs (1.1 kg)
Thigh Supports: 7.75"L x 5.75"W x 7.7"H (19.7 cm x 14.6 cm x 17.8 cm) • 1.4 lbs (0.6 kg)



Comfortable pads support chest, hips and thighs for long procedures.

- Multi-layered foam pads for optimal patient comfort.
- Easy-release locking clips for one-handed installation, adjustment and removal.
- Large hip supports available to accommodate larger patient anatomy.
- Intuitive design provides clear confirmation of locked or unlocked status at a glance.

DUAL RAIL COMPATIBLE† POSITION GUIDE

Curve™ Spine Frame System

KYRA9000K - Curve™ Complete Kit
KYRA9005K - Mini Curve™ Complete Kit

Pad set included.
Patient Weight Capacity: Standard - 600 lbs. (272 kg) • Mini - 500 lbs. (227 kg)
Device Weight: Standard - 24.4 lbs. (11 kg) • Mini - 22 lbs. (10 kg)

- Stable platform for spine procedures including laminectomies, decompressions, and disc surgeries.
- Intraoperative height adjustment for optimal surgical site exposure.

DUAL RAIL COMPATIBLE† POSITION GUIDE

Full product details, p. 54



WHAT'S INCLUDED



Connect Prone Support Set
Includes: Chest, 2x Hip, 2x Thigh (Left & Right)
KYRA9052K

OR



Connect Prone Support Set, Large
Includes: Chest, 2x Large Hip, 2x Thigh (Left & Right)
KYRA9053K

WHAT'S OPTIONAL



KYRA Comfort Covers
Disposable protection.
6 sets
KYRA9063 - Standard
KYRA9064 - Large

Why choose KYRA Connect Prone Supports?

Supports lock firmly in place—won't shift on rail



One-handed attach/detach with easy-release clips.



Easy-release locking clips allow for easy adjustment and maximum security.

† Compatible with Allen/Jackson-style operating tables with 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall rail sizes.
Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Find Your Fit

Model Comparison Guide






	Prone Support Set	Prone Support Set, Large	Curve	Mini Curve
COMPATIBILITY	Dual-rail table systems	Dual-rail table systems	Standard OR tables, dual-rail table systems and imaging tables.	Standard OR tables, dual-rail table systems and imaging tables.
WEIGHT CAPACITY	600 lbs. (272 kg)	600 lbs. (272 kg)	600 lbs. (272 kg)	500 lbs. (227 kg)
PATIENT SIZE	Optimal support for the average patient.	Larger hip supports provide support for larger patients.	For patients 5'7" and above.	For patients 5'9" and under.
USE CASE	Ideal for procedures that require flexible positioning and compatability with KYRA Connect and 'Jackson-style' tables.		Ideal for spine procedures requiring reverse lordosis including laminectomies, decompressions, fusions and discectomies.	

Arm Supports



Low-profile post enables use during proning without interference with imaging

Connect™ Swivel Armboard
KYRA2302

1" thick Clamshell Pad included. Sold individually.
Overall Dimensions: 16" x 7.5" x 4.6" (40.6 cm x 19.05 cm x 11.68 cm)
Armboard Weight: 7.0 lbs. (3.2 kg)
Patient Weight Capacity: 600 lbs. (272 kg)

- Designed for optimal positioning without C-arm interference.
- Lightweight, low-profile design.
- Clamshell Pad included for comfort and protection.
- Integrated SocketLok™ clamp for quick table attachment.

POSITION GUIDE

Specify Rail Type

US	KYRA2302-US	EU	KYRA2302-EU
UK	KYRA2302-UK	JP	KYRA2302-JP

WHAT'S INCLUDED



Connect Swivel Armboard
One armboard with clamshell pad

+



SocketLok™ Clamp
Permanently attached to armboard

WHAT'S OPTIONAL



Connect Rail Clamp, 6", p. 72
KYRA9525 (US)
KYRA9527 (EU)
KYRA9529 (JP)



No lost clamps—integrated rotating clamp allows multi-axis positioning



Silicone cover protects hardware from fluid and debris

A recommended set-up, using the
KYRA Connect™ Swivel Armboard



A **Connect™ Spine System**
KYRA9500

B **Prone Head Support System**
KYRA9071K

Prone Face Cushion, Adult
KYRA9080

C **Connect™ Head Plate**
KYRA9050

D **Connect™ Swivel Armboard**
KYRA2302

E **Connect™ Rail Clamp, 6"**
KYRA9525 (US)

F **Connect™ Prone Support Set**
KYRA9052K

Connect™ Comfort Covers
KYRA9063

G **Connect™ Table Mount Clamps**
KYRA9511 (US)

H **Connect™ Supine Top, Kit of 4** (not shown)
KYRA9056K

I **Connect™ Power Cord** (not shown)
KYRA9540 (US)

Get in touch to discover your ideal configuration
kyramedical.com

† Compatible with Allen/Jackson-style operating tables with 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall rail sizes.
Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Connect™ SPINE SYSTEM

Supine & Lateral Positioning

KYRA Connect™ Supine Top

Modular support pads for specialist spine systems.

FOR PATIENTS UP TO 600 LBS. (272 KG.)

The KYRA Connect Supine Top offers a simple way to adapt dual-rail tables for versatile positioning, creating a secure, even platform with layered cushioning designed to enhance patient comfort. Easy-release locking clips make setup, adjustments, and removal effortless. Available as a single top or in a convenient 4-piece kit.



One-handed attach/detach with easy-release clips.

Connect™ Supine Top

KYRA9056 - Single
KYRA9056K - Kit of 4

Pads included.
Patient Weight Capacity: 600 lbs. (272 kg)
Device Weight per section: 2.6 lbs. (1.2 kg)
Dimensions per section: 19"L x 12"W x 2.4"H (48.3cm x 30.5cm x 10.2cm)

- Modular design converts dual-rail tables into a stable, flat platform ideal for imaging.
- Base cutouts provide rail access for optimal accessory positioning.
- Easy-release locking clips for easy installation, adjustment and removal.
- Ultra-lightweight, user-friendly design.

DUAL RAIL COMPATIBLE† POSITION GUIDE



WHAT'S INCLUDED



Connect Supine Top
One supine top
KYRA9056

OR



Connect Supine Top, Kit of 4
Four supine tops
KYRA9056K

WHAT'S OPTIONAL



Connect Spine System Frame
KYRA9500



Easy-release locking clips for:

- Quick installation
- Seamless adjustments
- Simple removal

† Compatible with Allen/Jackson-style operating tables with 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall rail sizes.
Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Clamps & Accessories





DUAL RAIL COMPATIBLE†

Connect™ Rail Clamp, 6"

KYRA9525

Fits Accessory Rails: 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall
Patient Weight Capacity: 600 lbs. (272 kg)
Dimensions: 6"L x 4.5"W x 3.6"H (15.2 cm x 11.4 cm x 9.1 cm)

- ✓ Securely attaches to dual-rail systems to support attachment of standard OR table accessories (ex. armboards, clamps and retractors).

Specify Rail Type

US KYRA9525 EU KYRA9527 JP KYRA9529



DUAL RAIL COMPATIBLE†

Connect™ Rail Clamp, 12"

KYRA9526

Fits Accessory Rails: 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall
Patient Weight Capacity: 600 lbs. (272 kg)
Dimensions: 12"L x 4.5" W x 3.6"H (30.5 cm x 11.4 cm x 9.1 cm)

- ✓ Longer 12" rail clamp for greater flexibility in placement of standard OR table accessories (ex. armboards, clamps and retractors).

Specify Rail Type

US KYRA9526 EU KYRA9528 JP KYRA9530





DUAL RAIL COMPATIBLE†

Connect™ Rail Clamp, 12" Offset

KYRA9531 - Left
KYRA9532 - Right

Fits Accessory Rails: 1.25" (3.2 cm) wide and 1.5" (3.8 cm) tall
Patient Weight Capacity: 600 lbs. (272 kg)
Dimensions: 12"L x 3.9"W x 2.5"H (30.5 cm x 9.9 cm x 6.4 cm)

- ✓ Offset rail removes rail-space restrictions to allow optimal positioning of accessories (ex. retractors).

Specify Rail Type

US KYRA9531 (Left) KYRA9532 (Right)
EU KYRA9533 (Left) KYRA9534 (Right)
JP KYRA9535 (Left) KYRA9536 (Right)

† Compatible with Allen/Jackson-style operating tables with 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall rail sizes.
Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.



DUAL RAIL COMPATIBLE†

Connect™ Skull Clamp Adapter

KYRA9521

Sold individually. Universal compatibility.
Fits Accessory Rails: 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall
Patient Weight Capacity: 600 lbs. (272 kg)
Dimensions: 6"L x 4.5"W x 5.8"H (15.2 cm x 11.4 cm x 14.7 cm)
Product Weight: 4.4 lbs. (2 kg)

- ✓ Allows for secure attachment of Mayfield/ Doro-style skull clamps

Specify Rail Type

US KYRA9520 EU KYRA9521 JP KYRA9522

Connect™ Table Mount Clamps

KYRA9511

Includes left and right clamp.
Fits Accessory Rails: 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall
Dimensions: 8"L x 3.25"W x 2.5"H (20.3 cm x 8.3 cm x 6.4cm)

- ✓ Allows for attachment of Connect Spine System to your existing OR table.

Specify Rail Type

US KYRA9511 EU KYRA9512 JP KYRA9515

DUAL RAIL COMPATIBLE†





Connect™ Power Cord

KYRA9540

Dimensions: 10'5" long (3.2 m)

- ✓ Country-specific power cord for charging of Connect Spine System.

Specify Plug Type

US KYRA9540 EU KYRA9541
JP KYRA9542 UK KYRA9543



Connect™ Replacement Pendant Control

KYRA9544

- ✓ Replacement pendant for Connect Spine System.

ProneAlone™ Head Support

KYRA9082 - 8 per pack

Dimensions: 11.3"L x 9.7"W x 6.5" H (28.7 cm x 24.6 cm x 16.5 cm)
Product Weight: <1 lb. (<0.5kg)

- All-in-one standalone head support solution.
- Molded design for improved support, stability and alignment.
- Fully disposable including mirror for improved infection control.
- Mirror included for facial monitoring.

POSITION GUIDE



Disposable mirror included for patient observation



Foam Head Block

KYRA9022 - 6 per pack

Dimensions: 10.5"L x 9.75"W x 6"H (26.7 cm x 24.7 cm x 15.4 cm)
Product Weight: <1 lb. (<0.5kg)

- Provides support of the head and neck during prone positioning for 2 hours or less.
- Disposable for improved infection control.
- Slotted design to accommodate ventilation tubing.

POSITION GUIDE

CP Comfort Mask

KYRA9331 - 5 per pack

Dimensions: 11"L x12.75"W x 2"H (29.2 cm x 24.1cm x 5 cm)
Product Weight: <1 lb. (<0.5kg)

- Compatible with Allen® C-Flex™ Head Positioning System*.
- Contoured design helps relieve pressure around the eyes while supporting the face.
- Soft, cosmetic foam enhances patient comfort during extended procedures.

POSITION GUIDE



*Allen® and C-Flex™ are registered trademarks of Baxter.
Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Prone Face Cushion

KYRA9080 - Adult, 10 per pack
KYRA9083 - Child, 10 per pack

Dimensions: 12.7"L x 9.8"W x 6.1" H (32.3 cm x 24.9 cm x 15.5 cm)
Product Weight: 1.0 lbs. (0.5kg)

- For use with Prone Head Support System, p. 63.
- Unique molded design with safety rim; no sharp edges, improved protection during transfer.
- Evenly distributes pressure across the patient's face to reduce the risk of pressure-related injury.
- Adult and child cushions available to suit a wide range of patients.

POSITION GUIDE



WHAT'S INCLUDED



Prone Face Cushion, Adult
10 per pack
KYRA9080

OR



Prone Face Cushion, Child
10 per pack
KYRA9083

FOR USE WITH



Prone Head Support System, p. 63
KYRA9071K

Prone Face Cushion, Chinless

KYRA9081 - Adult, 10 per pack

Dimensions: 12.7"L x 9.8"W x 6.1" H (32.3 cm x 24.9 cm x 15.5 cm)
Product Weight: 1.0 lbs. (0.5kg)

- For use with Prone Head Support System, Chinless, p. 63.
- Unique molded design with safety rim; no sharp edges, improved protection during transfer.
- Chin-free profile preserves intubation access.
- Saves time by eliminating the need to remove and reinsert the airway during positioning.
- Improves patient safety by reducing contamination risk and avoiding unnecessary airway trauma.

POSITION GUIDE




WHAT'S INCLUDED



Prone Face Cushion, Chinless
10 per pack
KYRA9081

FOR USE WITH



Prone Head Support System, Chinless, p. 63
KYRA9072K



KYRA Comfort™ Covers

Disposable protection designed to improve infection control.

KYRA Comfort Covers deliver additional protection and support to both the patient and equipment, helping reduce skin shear and protecting hardware from excess fluids and wear and tear. The disposable design keeps room turnover simple and reduces contamination risk.



- Single-use to help reduce cross-contamination risk and simplify room turnover.
- Layer of supportive foam helps absorb excess moisture, dissipate heat and reduce risk of skin shear.
- Helps protects hardware from liquids and wear and tear.
- Sized to fit a variety of KYRA positioners.

Curve™ Comfort Cover Kit

KYRA9020

6 kits per pack

- For use with KYRA Curve, Bow and Wilson style frames.
- Each kit includes Curve Comfort Covers (2), Arm Cushions (2), Foam Head Block (1), Base Drape and Lateral Rod Pad (1).
- Curve Comfort Covers also sold separately below.





Connect™ Comfort Covers

KYRA9063 - For Prone Support Set
KYRA9064 - For Prone Support Set, Large

6 sets per pack

- For use with the Connect Prone Supports.
- Each set includes Chest Cover (1), Thigh Covers (2), Hip Covers (2).



Curve™ Comfort Covers

KYRA9020S - For Standard Curve
KYRA9024 - For Mini Curve

6 sets per pack

- For use with KYRA Curve, KYRA Mini Curve, Bow and Wilson style frames.
- Includes Curve Comfort Covers (2), Base Drape (1), Rod Support Cover (1).

Connect™ Comfort Covers, 3-Piece

KYRA9363 - For Prone Support Set
KYRA9364 - For Prone Support Set, Large

6 sets per pack

- For use with the Connect Prone Supports.
- Each set includes Chest Cover (1), Hip Covers (2).



Wingset Comfort Covers

KYRA9326

5 sets per pack

- For use with Wingset-style positioners.
- Each set includes 4 pad covers.



Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Curve™ Standard Cover Kit

KYRA9021 - 6 kits per pack

- Disposable support convenience kit for use with KYRA Curve, Bow and Wilson style frames.
- Single-use to help reduce cross-contamination risk and simplify room turnover.
- Helps protects hardware from liquids and wear and tear.
- Each kit includes Curve Standard Covers (2), Arm Cushions (2), Foam Head Block (1), Base Drape (1) and Lateral Rod Pad (1).





Curve™ Standard Cover

KYRA9021S - 6 per pack


- For use with KYRA Curve, Bow and Wilson-style frames.
- Single-use to help reduce cross-contamination risk and simplify room turnover.
- Helps protects hardware from liquids and wear and tear.

Arm Cushions

KYRA9023 - 6 pairs per pack

Dimensions: 23.9”L x 5”W x 3.8”H (60.7 cm x 12.7 cm x 9.6 cm)

- Supports the full length of the arm and hand, adjustable for different procedures.
- Anatomically contoured cutouts help maintain natural arm alignment and reduce pressure injury risk.
- Single-use pads help reduce cross-contamination risk.



Disposable Pad for Arthroscopic Leg Holder

KYRA8501 - 12 per pack

Dimensions: 37”L x 3.5”W x 1.25”H (94 cm x 8.9 cm x 3.2 cm)

- For use with the Arthroscopic Leg Holder, p. 48.
- Non-sterile disposable pad helps to protect from pressure and skin shear.
- Single-use pad for improved infection control.









Soft Disposable Armboard Strap

KYRA8310 - 12 per pack

Dimensions: 3”W x 32” L (7.6 cm x 81.3 cm)

- Disposable straps with hook-and-loop closure secure arm to armboard.
- Single-use for improved infection control.
- Minimal stretch ensures a firm hold while providing some conformity.

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Disposables

KYRA COMPLETE Commitment to Quality

At KYRA, we engineer our products for durability, performance, and long-term reliability.

Our warranty and service programs reflect the confidence we have in our materials, engineering, and manufacturing standards, and the commitment to our customers to protect your investment, minimize service disruption and keep your equipment operating at peak performance.

STANDARD WARRANTY

All KYRA products are covered by our standard manufacturer's warranty:

- **Hardware:** 2 years
- **Pads:** 1 year

Should an issue arise, our service team will respond promptly to assess and resolve it, ensuring minimal disruption to your clinical environment.

EXTENDED WARRANTY

KYRA stands behind its quality. Protect your equipment for maximum longevity with KYRAComplete™ Extended Warranty, our premium warranty and service program.

The KYRAComplete Extended Warranty includes:

- Up to 4 years of extended coverage
- A set of replacement pads (product dependent)
- Live troubleshooting support from our service team

KYRAComplete Extended Warranty provides enhanced protection and priority support, delivering predictable lifecycle management and added peace of mind.

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Maximum Uptime, Minimum Disruption

We recognize that reliability is essential in clinical environments. We strive to ensure equipment is easily repairable, pads replaceable, and support accessible for fast resolution if wear and tear or damage occurs—maximizing uptime and minimizing disruption to operating schedules.

On-site Repairs

Many common repairs can be completed by your Biomed department within your facility, eliminating the need for product returns and reducing operational disruption. With modular components and straightforward serviceability, KYRA solutions are built to keep you operational.

Easy-to-use Repair Kits

For added convenience, we offer dedicated repair kits for common service needs.

Each kit includes the required components and clear, step-by-step instructions, enabling your team to carry out repairs quickly and confidently. This approach maximizes uptime, avoids unnecessary shipping delays, and supports continuity of care.

Replacement Pads

Our pads are manufactured to high clinical and durability standards but wear and tear is inevitable over time.

We offer replacement pads across our full product range, designed for straightforward removal and installation. This allows you to refresh high-use surfaces while retaining your existing hardware, extending product lifespan and optimizing long-term value.

Replacement Parts & Servicing

Having an issue? Our service team is here to help.
Contact **service@kyramedical.com** and our technical specialists will:

- Provide guided troubleshooting
- Identify required replacement parts
- Support you through the repair process step by step

From initial inquiry to resolution, we are committed to responsive service and practical solutions that prioritize uptime and performance.

↑ Request a repair or service

Index

Order No.	Product Name	Page
KYRA0100	KYRAComplete Extended Warranty	80
KYRA1000	Purple Stirrup Cart	37
KYRA1010	Pediatric Stirrup Cart	37
KYRA1100	Curve Cart	55
KYRA1110	Mini Curve Cart	55
KYRA2000	Trigger Armboard, US	44
KYRA2005	Trigger Armboard, EU	44
KYRA2006	Trigger Armboard, UK	44
KYRA2007	Trigger Armboard, DY	44
KYRA2008	Trigger Armboard, JP	44
KYRA2009	Comfort Armboard Pad, 1"	45
KYRA2010	Comfort Armboard Pad, 2"	45
KYRA2015	Comfort Armboard Pad, 3"	45
KYRA2050	Clamshell Armboard Pad, 1"	45
KYRA2051	Clamshell Armboard Pad, 2"	45
KYRA2052	Clamshell Armboard Pad, 3"	45
KYRA2060	Bariatric Armboard	43
KYRA2300	Swivel Armboard, Large	46
KYRA2304	Swivel Armboard, Small	46
KYRA2305	Side Rail Extender, US	41
KYRA2306	Side Rail Extender, EU	41
KYRA2320	Table Width Extender	42
KYRA2321	Table Width Extender Pad	42
KYRA2330	Arm Protection Shield	43
KYRA2340	Secondary Shoulder Supports	51
KYRA2401	Patient Transfer Board, Long/Wide	49
KYRA2403	Patient Transfer Board, Short/Narrow	49
KYRA2404	Patient Transfer Board, Short/Wide	49
KYRA2501	Flexible Anesthesia Screen	42
KYRA2601	Patient Safety Strap, Medium	51
KYRA2602	Nissen Straps, Bariatric	50
KYRA2603	Nissen Straps, Standard	50
KYRA3000	Comfort 350 Stirrup, Hook & Loop	20-21
KYRA3100	Comfort 350 Stirrup, Secure-Lok	20-21
KYRA3200	Clean 350 Stirrup	28-29
KYRA3800	Classic 350 Stirrup	10-12
KYRA4100	Secure Release Clamp, US	40
KYRA4200	Secure Release Clamp, EU	40
KYRA4300	Secure Release Clamp, UK	40
KYRA4400	Secure Release Clamp, DY	40
KYRA4500	Secure Release Clamp, JP	40
KYRA4650	BladeLok Clamp, US	40
KYRA4660	BladeLok Clamp, EU	40
KYRA4670	BladeLok Clamp, UK	40
KYRA4690	BladeLok Clamp, JP	40
KYRA4700	Simple Blade Clamp, US	40
KYRA4710	Simple Blade Clamp, EU/UK	40
KYRA4730	Simple Blade Clamp, DY	40
KYRA4740	Simple Blade Clamp, JP	40
KYRA4850	SocketLok Clamp, US	41
KYRA4860	SocketLok Clamp, EU	41
KYRA4870	SocketLok Clamp, UK	41
KYRA4890	SocketLok Clamp, JP	41
KYRA5000	Comfort 500 Stirrup, Hook & Loop	22-23
KYRA5100	Comfort 500 Stirrup, Secure-Lok	22-23
KYRA5200	Clean 500 Stirrup	30-31
KYRA5800	Classic 500 Stirrup	12-13
KYRA5850	Classic 500 Plus Stirrup	12-13
KYRA6200-P	Comfort 500 Repl. Pads, Hook & Loop	19
KYRA6210-P	Comfort 350 Repl. Pads, Hook & Loop	17
KYRA6250-P	Comfort 500 Repl. Pads, Secure-Lok	19
KYRA6260-P	Comfort 350 Repl. Pads, Secure-Lok	17
KYRA6270-P	Clean 500 Repl. Pads	27
KYRA6290-P	Clean 350 Repl. Pads	25
KYRA6320-P	Comfort 800 Repl. Pads	21
KYRA6350-P	Classic 350 Repl. Pads	11
KYRA6360-P	Classic 500 Plus Repl. Pads	13
KYRA6370-P	Classic 500 Non-Fin Repl. Pad	15
KYRA6380-P	Classic 800 Plus Repl. Pad	17
KYRA6470-P	Classic 500 Repl. Pads	13
KYRA7100	KYRA 400 Stirrup	32
KYRA7300	Knee Crutch, Goepel	38
KYRA7305	Knee Crutch, Lift Assist	33
KYRA7800	Classic 500 Non-Fin Stirrup	14-15
KYRA7810	KYRA 250 Stirrup	32
KYRA7860	Simple Stirrups	38
KYRA8100	Comfort 800 Stirrup	24-25
KYRA8310	Soft Disposable Armboard Straps	79
KYRA8400	Urology Attachment Band, Universal	36
KYRA8410	Urology Drainage Bags, Non-Sterile	36
KYRA8411	Urology Drainage Bags, Sterile	36
KYRA8451	Clean Catcher Drainage System, Sterile	34
KYRA8452	Clean Catcher Drainage System, Non-Sterile	34
KYRA8460	Clean Catcher Stirrup Mounts, Round	35
KYRA8461	Clean Catcher Stirrup Mounts, Flat	35
KYRA8462	Clean Catcher Stirrup Mounts, Allen	35
KYRA8501	Disposable Pad for Arthroscopic Leg Holder	48
KYRA8850	Classic 800 Plus Stirrup	16-17
KYRA9000	Curve Spine Frame	54-55
KYRA9000K	Curve Complete Kit	54-55
KYRA9005	Mini Curve Spine Frame	54-55
KYRA9005K	Mini Curve Complete Kit	54-55
KYRA9010	Curve Crank Handle	55
KYRA9020	Curve Comfort Covers Kit	77
KYRA9020S	Curve Comfort Covers	77
KYRA9021	Curve Standard Covers Kit	78
KYRA9021S	Curve Standard Covers	78
KYRA9022	Foam Head Block	74
KYRA9023	Arm Cushions	78
KYRA9024	Mini Curve Comfort Covers	77
KYRA9030	Curve Strap Buckle	55
KYRA9050	Connect Prone Head Plate	64
KYRA9052K	Connect Prone Support Set	66
KYRA9053K	Connect Prone Support Set, Large	66
KYRA9056	Connect Supine Top	70-71
KYRA9056K	Connect Supine Top, Kit of 4	70-71
KYRA9063	Connect Comfort Covers	76
KYRA9064	Connect Comfort Covers, Large	76
KYRA9067	Connect Pivoting Mirror	64
KYRA9071	Prone Helmet	62-63
KYRA9071K	Prone Head Support System	62-63
KYRA9072	Prone Helmet, Chinless	62-63
KYRA9072K	Prone Head Support System, Chinless	62-63
KYRA9080	Prone Face Cushion, Adult	75
KYRA9081	Prone Face Cushion, Chinless, Adult	75
KYRA9082	ProneAlone Head Support	74
KYRA9083	Prone Face Cushion, Child	74
KYRA9100K	Lateral Positioning Set	52-53
KYRA9110	Lateral Positioner, Round	53
KYRA9120	Lateral Positioner, Rectangle	53
KYRA9200	Arthroscopic Leg Holder	48
KYRA9205	Well Leg Holder	39
KYRA9326	Wingset Comfort Covers	77
KYRA9331	CP Comfort Mask Covers	74
KYRA9363	Connect Comfort Covers, 3-Piece	76
KYRA9364	Connect Comfort Covers, Large, 3-Piece	76
KYRA9500	Connect Spine System Frame	56-60
KYRA9511	Connect Table Mount Clamps, US, L/R	73
KYRA9512	Connect Table Mount Clamps, EU, L/R	73
KYRA9515	Connect Table Mount Clamps, JP, L/R	73
KYRA9521	Connect Skull Clamp Adapter	73
KYRA9525	Connect Rail Clamp, 6", US	72
KYRA9526	Connect Rail Clamp, 12", US	72
KYRA9527	Connect Rail Clamp, 6", EU	72
KYRA9528	Connect Rail Clamp, 12", EU	72
KYRA9529	Connect Rail Clamp, 6", JP	72
KYRA9530	Connect Rail Clamp, 12", JP	72
KYRA9531	Connect Rail Clamp, 12" Offset, US, Left	72
KYRA9532	Connect Rail Clamp, 12" Offset, US, Right	72
KYRA9533	Connect Rail Clamp, 12" Offset, EU, Left	72
KYRA9534	Connect Rail Clamp, 12" Offset, EU, Right	72
KYRA9535	Connect Rail Clamp, 12" Offset, JP, Left	72
KYRA9536	Connect Rail Clamp, 12" Offset, JP, Right	72
KYRA9540	Connect Power Cord, US	73
KYRA9541	Connect Power Cord, EU	73
KYRA9542	Connect Power Cord, JP	73
KYRA9543	Connect Power Cord, UK	73
KYRA9544	Connect Replacement Pendant Control	73

KYRA MEDICAL

MADE IN THE USA

REINVENTING THE CLASSIC

KYRA® Classic™ Stirrups

Combining patented KYRA technology with classic stirrup design—the Classic stirrup features an improved easy-squeeze handle and traditional straight rod, and can treat your tallest to shortest patients.



Fabric-free silicone straps
help improve infection control.

Easy-squeeze Handle
Reduced squeeze angle and grip force suits smaller hands.

Unique Design
Patented design supports safe positioning of the tallest to shortest patients.

Learn more about the new **Classic™ Stirrups**
p. 8

Shown with KYRA Clean™ Catcher Drainage System, p. 34

KYRA MEDICAL

Your trusted partners in patient positioning

KYRA Medical, 102 Otis Street, Northborough, MA 01532 USA
kyramedical.com | +1-508-936-3550

©2026 KYRA Medical, Inc. All Rights Reserved. KYRA0010-A12